English

A. DEVICE DESCRIPTION

The FibuLock® Nail is designed to provide fixation of fibula fractures while they heal.

B. INDICATIONS

The FibuLock System is intended for use in the fixation of fibula fractures and osteotomies.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone/soft tissue.
- Patients with poor prognosis for healing.
- Blood supply limitations and previous infections, which may retard healing.
- Nickel sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Active or latent infection.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Potential adverse events related to any surgery: Reactions to anesthesia, the anesthetic or other medications, bleeding, infection (both deep and superficial), ileus, blood vessel damage, nerve or soft tissue damage, atelectasis, pneumonia, hematoma, seroma, wound dehiscence or incisional hernia, urologic problems, embolism, anemia, colitis, thrombophlebitis, heart attack, stroke or death.
- Potential adverse events related to the device: Bending, cracking, or fracture of the implant, failure to achieve fracture healing, or loosening of the implant. Limb shortening or loss of anatomic position with nonunion, malunion with rotation or angulation. Irritational injury of soft tissues, including impingement syndrome, or tissue reactions which include macrophage and foreign body reactions adjacent to implants, Metal sensitivity reactions and/or allergic reactions to foreign materials, pain, discomfort, or abnormal sensations due to the presence of the implant. Necrosis of bone or bone resorption. Necrosis of the tissue or inadequate healing. Generation of material debris possibly resulting in bodily response. Additional surgery may be necessary for implant removal, repositioning or replacement.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- Do not re-sterilize this device.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Bone fixation devices are neither intended to carry the full load of the patient, nor intended to carry a significant portion of the load for extended periods of time. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- The correct selection and placement of the implant is extremely important. The appropriate type and size should be selected for the patient. Failure to use the correct implant size or improper positioning may result in loosening, bending, cracking, or fracture of the device, or bone or both.

- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- The patient must be cautioned about the use, limitations, and possible adverse effects of this implant including the possibility of the device failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight and/or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and electromagnetc simulations demonstrated that the FibuLock Nail Ankle Fracture System is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the FibuLock Nail Ankle Fracture System is expected to produce a maximum temperature rise less than 6 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the FibuLock Nail Ankle Fracture System extends approximately 40 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Implants are single-use only and should never be re-used.
- Use care in handling and storage of implant components. Cutting, sharply bending or scratching the surface can significantly reduce the strength and fatigue resistance of the implant system.
- The surgeon should be thoroughly familiar with the devices, instruments and surgical technique before surgery. A detailed surgical technique is available for surgeon reference. Proper surgical procedures and techniques are the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of the procedures based on his or her personal medical training and experience. During the post-operative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed about post-surgical activity requirements.
- Implants should not be forced, twisted, bent, or pried during insertion as this may damage, break, or bend the implant and render it unusable. Damaged implants should not be implanted.
- The use of locking screws is necessary for fracture stability. Refer to the surgical technique or product listing for the available screw sizes.
- The surgeon should take care with patient selection when treating patients with preexisting conditions that may affect the success of the surgical procedure.
- The patient should be advised that an additional procedure for the removal of the implant is sometimes necessary. Providing postoperative instructions to patients and appropriate nursing care is critical. Early weight and/or load bearing substantially increases implant loading and increases the risk of loosening, bending or breaking the device. Patients who are obese and/or noncompliant, as well as patients who could be pre-disposed to delayed union or non-union, must have auxiliary support. Patients should be cautioned against unassisted weight-bearing activity without physician direction or medical release. Postoperative care and physical therapy should be structured to prevent loading of the operative extremity until directed by the physician.
- Even after full healing, the patient should be cautioned that re-fracture is possible and that a revision surgery may be necessary.
- Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The FibuLock Nail is manufactured from medical grade stainless steel.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. REVISION SURGERY / IMPLANT REMOVAL

Metallic implants can loosen, fracture, corrode, migrate, cause pain, or stress shield bone even after a fracture is healed. If any of these complications occur, the surgeon must make the final decision on implant removal.

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Arthrex, Inc.
Arthrex GmbH