Symbols Used on Labeling
 Auf der Verpackung verwendete Symbole
 Symboles utilisés sur l'étiquetage
 Simboli usati sull'etichetta
 Símbolos utilizados en las etiquetas
 Símbolos utilizados na embalagem

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English DFU-0197r1

A. DEVICE DESCRIPTION
 The Single Use Harvester for Osteochondral Repair is a size specific cylindrical Harvester with a core driver and Sizer/Tamp rod.

website also provides detailed surgical technique information and demonstrations.

B. INDICATIONS
 The Harvester is used for the removal of an osteochondral cylindrical block, which is replaced by a repair implant of the corresponding size.

C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.
3. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
4. Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Increased effusion or morbidity.

E. WARNINGS

1. Arthrex recommends manual-skills training on cadaveric or animal specimens with associated instrumentation prior to the procedure. The surgical technique manual should be reviewed prior to the procedure.
2. All potential harvest site locations should be researched and understood prior to the procedure.
3. The length and diameter of the harvested osteochondral tissue MUST precisely match the length and diameter of the implant which will be used to fill the defect.
4. This device must never be reused.
5. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

F. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.

H. STERILIZATION
 This device is supplied provided sterile. Refer to the package label for the sterilization method.

I. STORAGE CONDITIONS
 Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

J. INFORMATION
 Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques. Or, contact your Arthrex representative for an onsite demonstration. The Arthrex

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