English

A. DESCRIPTION

The **FiberTape® Cerclage** and **TigerTape™ Cerclage** sutures are assembled from **FiberTape®** and **TigerTape™** sutures, in a loop configuration. The looped tip of the suture is tied as a hitch over a suture sheath securing a shuttling polyblend suture.

FiberTape and TigerTape are non-absorbable, flat braided sutures composed of UHMWPE and polyester yarns over a core of FiberWire® or TigerWire® sutures (each made of a polyblend of UHMWPE and polyester). Additional materials may include cyanoacrylate and/or nylon.

These sutures meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements).

These sutures are supplied sterile, in pre-cut lengths, in various single and multi-loop configurations, and in some cases, with various swaged needles, with stiffened ends and with a pre-tied hitch. The sutures are available non-dyed, dyed, or fully or partially striped. Dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon.

B. INDICATIONS

FiberTape Cerclage, TigerTape Cerclage, FiberTape and TigerTape sutures are intended for use in soft tissue approximation and or ligation. These sutures may be incorporated, as components, into surgeries where constructs including those with allograft or autograft tissues are used for repair.

When used as bone fixation cerclage the sutures are intended for

- Trochanteric reattachment after trochanteric osteotomy following total hip arthroplasty.
- Sternotomy indications including the “rewiring” of osteomized sternums.
- Trauma surgery indications including olecranon, ankle, patella and some shoulder fracture rewiring.
- Repair of long bone fractures due to trauma or reconstruction.

C. ACTIONS

These sutures when tested per ISO 10993, *Biological Evaluation of Medical Devices-Part 10: Tests for Irritations and Sensitization*, had no reactions of allergic or sensitive nature. The dyed suture and coating are pharmacologically inactive.

These sutures are not absorbed, but may become encapsulated in the surrounding connective tissues. The sutures are not known to have significant change in tensile strength *in vivo*.

D. CONTRAINDICATIONS

- These sutures are not for use in cardiac indications.
- Spinal fixation applications.
- Osteoporosis or insufficient bone quality.
- Patients that are skeletally immature.
- Cases where there is poor quality or osteoporotic bone at the repair.
- Compromised vascularity that would inhibit adequate blood supply to the fracture or operative site.

E. WARNINGS

- Do not re-sterilize.
- Once open, discard unused suture.
- Do not expose suture to heat.
- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Arthrex sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.
- As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formations. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- The safety and effectiveness of these devices for use as an artificial ligament or tendon has not been established; hence, the device should not be used for those applications.

F. MRI SAFETY INFORMATION

1. This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.

G. PRECAUTIONS

- In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- Assure that all knots have been secured using ac- cepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particu- larly appropriate when knotting monofilaments. Care should be taken to prevent damage to surrounding tissue or user puncture due to improper handling of the needlepoint.
- Do not grasp the needle at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bend- ing and breaking. Discard used needles in “sharps” containers.
- Avoid wrapping the suture over sharp metal or bone graft surfaces.

H. ADVERSE REACTIONS

- Adverse reactions have not been noted with the Arthrex FiberTape product in animal testing.
- Common non-absorbable suture reactions may include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site.
- Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.
- Loosening.
- Need for revision surgery.

I. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/ symbolsglossary](http://www.arthrex.com/symbolsglossary).

J. STERILIZATION

The device is provided sterile. Refer to the package label for the sterilization method.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79, “Com- prehensive Guide to Steam Sterilization and Sterility As- surance in Health Care Facilities” for specific information.

K. STORAGE CONDITIONS

This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

Surgeons are advised to review the product specific surgi- cal technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstra- tions. Or, contact your Arthrex representative for an onsite demonstration.

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