English

A. DEVICE DESCRIPTION

These devices consist of metal FASTak® and Corkscrew® Suture anchors.

The FASTak and Corkscrew suture anchors are threaded implants with eyelets.

B. INDICATIONS

The FASTak and the Corkscrew suture anchors are intended to be used for suture (soft tissue) fixation to bone in the foot, ankle, knee, hand, wrist, elbow, shoulder and hip (devices with FiberWire® only).

Surgeons must apply their professional judgment when determining the appropriate suture-anchor size based on the specific indication, preferred surgical technique, and patient history.

These products are intended for the following indications:

Shoulder: Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulolabral Reconstruction.

Foot/Ankle: Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Mid-foot Reconstruction, Metatarsal Ligament Repair and (Corkscrew suture anchor only) Bunionectomy.

Knee: Medial Collateral Ligament Repair, Lateral Colateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis and (Cork-screw suture anchor only) Anterior Cruciate Ligament Repair.

Hand/Wrist: Scapholunate Ligament Reconstruction, Repair/Reconstruction of Collateral Ligaments.

Elbow: Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction and (Corkscrew suture anchor only) Biceps Tendon Reattachment, Lateral Epicondylitis Repair

Hip: Capsular repair, acetabular labral repair.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or will- ingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immat- ure. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Shoulder dislocation/subluxation.

E. WARNINGS

- An internal fixation device must never be reused.
- All metallic implant devices used for this surgical proce- dure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important consid- erations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- This device has not been evaluated for safety and com- patibility in the magnetic resonance (MR) environment. This device has not been tested for heating or migration in the MR environment. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.
- Appropriate selection of size and type of suture prior to use of these devices is advised. Arthrex recommends use only with FiberWire®, TigerWire®, FiberTape®, and TigerTape™ sutures **or** 100% Ultra High Molecular Weight Polyethylene (UHMWPE).

F. MRI SAFETY INFORMATION

I. MR CONDITIONAL

Non-clinical testing and in-vivo electromagnetic simula- tions demonstrated that the Economy Anchors are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla. only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body aver- aged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Economy Anchors are expected to produce a maximum temperature rise of up to 1.8°C after 15-minutes of continuous scanning.

- a. Artifact Information**

In non-clinical testing, the image artifact caused by the Economy Anchors can extend up to approximately 17 mm from this implant when imaged using a gradient ecbo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolglossary.

I. STERILIZATION

This device is provided non-sterile and can be sterilized. It must be adequately cleaned, then sterilized using one of the following sterilization parameters.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
Gravity-Displacement Steam Sterilization Cycle	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C (270°F)	15 Minutes	15 to 30 Minutes
	135°C (275°F)	10 Minutes	30 Minutes
Pre-vacuum Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
Gravity-Displacement Steam Sterilization Cycle	132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
	Pre-vacuum Cycle	132°C - 135°C (270°F - 275°F)	4 Minutes

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities” for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions. Cooling – The device must be adequately cooled, after being removed from the sterilizer.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The device is made of titanium.

K. STORAGE CONDITIONS

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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