## English

### A. DEVICE DESCRIPTION

The PIP (Proximal interphalangeal) Dart is an implant with a barbed design and ranges in size from 20 to 35 mm in length. The PIP Dart is inserted between the proximal and middle phalanges, so the barbs fixate on the phalangeal canal of the toe.

The RetroFusion™ screw is an implant with a dual threaded design and ranges in size from 20 to 32mm in length. The RetroFusion screw is inserted between the proximal and middle phalanges, so the opposing threads fixate on the phalangeal canal of the toe and compress the joint.

The TRIM-IT Drill Pin® is a disposable kit containing a Bioabsorbable implant, K-Wire, Tamp and Guide Sleeve. The implant has a diameter of 1.5 mm or 2 mm and length of 100 mm. It is pointed at both ends.

The TRIM-IT Spin Pin™ is a disposable kit containing an implant and ruler. The implant has a diameter of 2 mm with a length of 102 mm. One end is made of Bioabsorbable material, while the other end has a metal tip. The metal tip has graduated lines for measuring. Both ends are pointed.

The Chondral Dart™ facilitates fixation and compression of osteochondral flap tears up to 2 cm (0.79 in.) in diameter.

### B. INDICATIONS

The PIP Darts and RetroFusion screws are intended to fix small bony or apical chondral fragments in the foot, ankle, upper extremities, hand, and wrist, where such fragments are not under tension or load-bearing. These devices are used in cases of osteochondritis dissecans and osteochondral fragments, fixation of fractures, 1<sup>st</sup> metatarsal (bunionectomy osteotomies), cuneiform bones, inherently stable osteotomies, and fusions of the phalanges, metatarsals, metacarpals, carpal bones, tarsal bones, ankle, and wrist. The devices can be used for inherently stable intramedullary stabilization of joint arthroplasty (resection) or fusion for the treatment of digital deformities of the foot or hand. This device is also used in inherently stable long bone fractures such as the femur, fibula, tibia, radius and ulna, including the diaphyseal, epiphyseal, and metaphyseal areas.

The Drill Pin and the Spin Pin are intended for fixation of small bony or apical chondral fragments in the Foot and Ankle where such fragments are not under tension or load-bearing. These pins are used in cases of osteochondritis dissecans and osteochondral fragments, fixation of fractures, 1<sup>st</sup> metatarsal (bunionectomy osteotomies), cuneiform bones, inherently stable osteotomies, and fusions of the phalanges, metatarsals, tarsal bones, and ankle. The Drill Pin and the Spin Pin can be used for inherently stable intramedullary stabilization of joint arthroplasty (resection) or fusion for the treatment of digital deformities of the foot.

The Chondral Dart is intended for the use in the fixation of small bone fragments, such as apical fragments, osteochondral fragments, and cancellous fragments in the knee.

### C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may tend to retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

### D. ADVERSE EFFECTS

- Infection, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

### E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be reused.
- Over-flexing of the device may cause breakage or cracking.
- Postoperatively, until healing is complete the fixation provided by this device should be protected. The post-operative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the implant, otherwise loosening, fracture, or migration of the device may result.

- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Metal only:** Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- The 1.5mm Drill Pin can be used with the hammer toe indication only.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
- Follow your institutions policy for safe disposal of all needles and other sharps or medical waste.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

### F. MRI SAFETY INFORMATION

#### 1. MR Conditional

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the metal (titanium and stainless steel) PIP Darts, RetroFusion Screws and TRIM-IT Devices are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the PIP Darts, RetroFusion Screws and TRIM-IT Devices are expected to produce a maximum temperature rise of up to 14.9 °C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the PIP Darts, RetroFusion Screws and TRIM-IT Devices can extend up to approximately 120 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

- MR Safe**
- The PIP Darts, RetroFusion Screws, TRIM-IT Devices, and Chondral Darts manufactured from only poly (L-Lactide) and/or polyetheretherketone (PEEK) are MR safe.*

### G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

### H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbolsglossary**.

### I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

- Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

### J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The device is made of Polyetheretherketone (PEEK), Titanium or Poly L-lactide (PLLA).

The TRIM-IT Spin Pin has a metal trocar tip attached to one end.

### K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-Bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

### L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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## PIP Darts, RetroFusion™ Screws, TRIM-IT™ Devices, and Chondral Darts

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This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the “Arthrex U.S. Product Warranty” section of the Arthrex, Inc. website, found at www.arthrex.com whose provisions are incorporated herein by reference.

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