<p><i>Symbols Used on Labeling</i> Auf der Verpackung verwendete Symbole Symboles utilisés sur l’etiquette Simboli usati sull’etichetta Símbolos utilizados en las etiquetas Símbolos utilizados na embalagem</p>

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English DFU-0214r1

A. INDICATIONS

MultIGun is a handheld non-invasive reusable instrument intended for use with Arthrex Quickset™ syringes.

B. WARNINGS

- Prior to use, inspect and verify that the device works properly (see section Q).
- Do not use to manipulate tissue or for any use other than indicated.
- Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
- The MultiGun should only be used by trained personnel. There is a risk of hand injury or damage to the protective glove in the case of improper handling. Avoid catching fingers when attaching the dual-chambered Quickset™ syringe.

C. LIMITATIONS ON REPROCESSING

Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.

D. VALIDATION

The recommended cleaning, disinfection, and sterilization methods in this DFU have been validated in compliance with federal and international guidelines/standards. Cleaning, disinfecting, and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility/end user to perform the appropriate validation testing.

- The following instructions were developed using the guidance given in the following standards:
- ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities"
 - ISO 17664: Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices

E. CONTAINMENT AND TRANSPORTATION

It is recommended that instruments are reprocessed as soon as is reasonably practical following use.

F. PREPARATION FOR CLEANING

When properly performed, cleaning, disinfection and/or sterilization do not compromise the use and mechanical performance of these instruments. These instruments are used with or on patients who may harbor both recognized and unrecognized infections. To prevent the spread of infection, all reusable instruments must be thoroughly cleaned, disinfected, and sterilized after use on each patient.

- Separate the Dispenser Piston from the MultiGun prior to cleaning.
- Remove dried-on soil from devices, especially in areas such as joints and crevices, prior to washing.

G. INSPECTION AND MAINTENANCE

- Arthrex non-sterile instruments are precision medical instruments and must be used and handled with care.
- Inspect the instruments for damage prior to use and at all stages of handling thereafter.
- Devices with cutting functions or sharp points become dull with continuous use. This condition does not indicate a device defect. This condition indicates normal wear. Dull devices may require replacement if they no longer perform as designed. Inspection prior to use should include verifying the cutting ability and sharpness of these edges.
- If damage is detected, do not use the device prior to consulting the manufacturer for guidance.
- Dry instruments thoroughly and lubricate all moving parts with a water-soluble instrument lubricant prior to sterilization.
- Check instruments for visible soil. Repeat cleaning if soil is visible and reinspect.

H. MANUAL CLEANING

- Immediate rinsing and cleaning after use with an enzymatic or alkaline cleaning detergent will effectively remove and prevent drying of adherent blood, mucus, etc. Cleaning solutions can include, but are not limited to: ENZOL® enzymatic, Neodisher® Mediclean Forte, and Thermosep® alkali clean. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cycolac™.**
- Scrub instrument with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the instrument.
- Rinse the instrument thoroughly with water following the cleaning process.
- Check instruments for visible soil. Repeat cleaning if soil is visible and reinspect.

I. ULTRASONIC CLEANING

- The instrument should be placed in an ultrasonic cleaning unit for a minimum of 20 minutes and processed according to the ultrasonic unit's directions.
- The instrument should be rinsed thoroughly with water following the ultrasonic process.
- Check instruments for visible soil. Repeat cleaning if soil is visible and reinspect.

J. AUTOMATIC WASHING

- Load the instruments in the washer such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain (fingers should be open and cannulations/holes positioned to drain).
- Run the automatic wash cycle – minimum cycle parameters:
 - 2 minute cold prewash at 68 ± 9°F (20 ± 5°C).
 - 3 minute cleaning wash (enzymatic or alkaline agent) at 140 ± 9°F (60 ± 5°C).
 - 15 second rinse at 140 ± 9°F (60 ± 5°C).
 - 1 minute thermal rinse at 176 ± 9°F (80 ± 5°C).
 - 6 minute drying phase at high temperature.
- Automatic wash cleaning solutions can include, but are not limited to: ENZOL® enzymatic, Neodisher® Mediclean Forte, and Thermosep® alkali clean. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cycolac™.**
- Check instruments for visible soil. Repeat cleaning if soil is visible and reinspect.

K. MANUAL DISINFECTION

- Instruments should be cleaned before disinfection, as blood albumen will impair the bactericide effectiveness of the solution. Immense instruments in disinfection solutions for a minimum of 20 minutes.
- Suitable disinfection solutions can include, but are not limited to: CIDEX® WAI-CIDE®-01, Gigasept®, kohrsolin® and equivalent products). Use the supplier's instructions for preparing the solution. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cycolac™.**
- After disinfection, the instruments should be rinsed with distilled water or preferably demineralized sterile water.
- Dry instruments thoroughly and lubricate all moving parts with a water-soluble, medical-instrument lubricant prior to sterilization.

L. STERILIZATION

This device is provided non-sterile.
Certain Arthrex devices that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use.
Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.

Cooling – The instrument must be adequately cooled, after being removed from the sterilizer. It should not be touched during the cooling process. Do not place the instrument on a cold surface or immerse in a cold fluid.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:

	Exposure Temperature	Exposure Time	Drying
			Time
Gravity-Displacement Steam Sterilization Cycle	132°C (270°F)	15 Minutes	15 to 30 Minutes
	135°C (275°F)	10 Minutes	30 Minutes
Pre-vacuum Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:

	Exposure Temperature	Exposure Time	Drying
			Time
Gravity-Displacement Steam Sterilization Cycle	132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes

M. PACKAGING

Simply a standard packaging material may be used. Ensure that the pack is large enough to contain the instrument without stressing the seals.

Sets: Instruments may be loaded into dedicated instrument trays or general-purpose sterilization trays. Ensure that cuttings edges are protected and do not exceed 8.5 kg/18.7 lb per tray. Wrap the trays using the appropriate method.

N. STORAGE

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

O. SPECIAL PRECAUTION - TRANSMISSIBLE SPONGIFORM ENCEPHALOPATHY AGENTS

It is outside the scope of this document to describe in detail the precautions that should be taken for Transmissible Spongiform Encephalopathy Agents. The agents for transmission of Creutzfeldt-Jakob disease are believed to be resistant to normal processes of disinfection and sterilization and therefore the normal processing methods of decontamination and sterilization as described above may not be appropriate where CJD transmission is a risk.

In general, the issues that come into contact with orthopedic surgical instruments are those of low TSE infectivity. However, particular precautions should be taken when handling instruments that have been used on known, suspected, or at-risk patients.

P. TECHNICAL SPECIFICATIONS

- Weight: approx. 730 g
- Power transmission ratio: 6:1
- Stroke length Dispenser Piston: > 3.5 mm
- Break load on Dispenser Piston: < 3000 N

Q. INSPECTION

- Before use the device should be inspected for signs of wear or damage.
- In particular look for any of the following problems:
- Damage to the Syringe Adapter
 - Damage to the Trigger
 - Damage to the Catch
 - Damage to the Dispenser Piston
 - Obvious damage to the unit
 - Traces of blood, debris, stains, and bone substitute

R. DIRECTIONS FOR USE

Element Number	Element Name
1	Syringe Adapter
2	Trigger
3	Catch
4	Dispenser Piston

- Insert the prepared dual-chambered Quickset™ syringe (please refer to the current version of the Arthrex Quickset™ guide) into the Syringe Adaptor of the MultiGun with the powder chamber facing downwards.
- Unlock the Catch by counter-clockwise rotation to .
- Insert the Dispenser Piston through the Catch and push the Dispenser Piston forward until it gets in contact with the dual-chambered Quickset™ syringe.
- Lock the Catch by clockwise rotation to .
- Ensure that the collar of the dual-chambered Quickset™ syringe is turned to "inject".
- If needed, connect the 7G cannula to the luer tip of the Quickset™ syringe.
- Compress the Trigger of the MultiGun in order to push the Dispenser Piston forward and expel any excess air. Dispense Arthrex Quickset™ resorbable bone substitute by compressing the Trigger. Assure the application of Arthrex Quickset™ within 4 min after initial transfer of the liquid phase into the powder chamber during the Quickset™ syringe preparation (please refer to the current version of the Arthrex Quickset™ guide).
- On completion of application, unlock the Catch of the MultiGun by counter-clockwise rotation to .
- Pull the Dispenser Piston back and remove the dual-chambered Quickset™ syringe. Do not compress the Trigger of the MultiGun while retracting the Dispenser Piston.

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