English

A. DEVICE DESCRIPTION

The Angel® PRP (Platelet Rich Plasma) Kit contains blood draw components inclusive of an infusion set, tourniquet, alcohol pad, gauze, Band-Aid, 60mL syringes, in addition to the Angel cPRP (concentrated Platelet Rich Plasma) Processing Set.

B. INDICATIONS

In the USA, The Arthrex Angel System Kits are to be used intraoperatively at the point of care for the safe and rapid preparation of autologous platelet poor plasma and platelet concentrate (platelet rich plasma) from a small sample of peripheral blood or a small sample of a mixture of peripheral blood and bone marrow. The platelet poor plasma and platelet rich plasma are mixed with autograft and/or allograft bone prior to application to a bony defect for improving handling characteristics.

Disclaimer – Platelet Rich Plasma prepared from a mixture of whole blood and bone marrow may contain higher levels of plasma free hemoglobin than Platelet Rich Plasma prepared from whole blood.

Outside the US, the Angel PRP Kit is intended to facilitate blood draws for use with the Angel cPRP Processing Set in the production of PRP.

The Angel cPRP Processing Set is intended for use with the Angel Concentrated Platelet Rich Plasma (cPRP) System to separate and collect an autologous plasma fraction rich in platelets and red blood cells from the patient's whole blood or a small mixture of blood and bone marrow perioperative to a surgical procedure, for treatment of osteoarthritis of the knee or rotator cuff repair.

C. CONTRAINDICATIONS

1. Insufficient quantity or quality of bone.
2. Blood supply limitations and previous infections, which may retard healing.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
4. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
5. Do not use for surgeries other than those indicated.
6. Any active infection or blood supply limitations.
7. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
8. Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Foreign body reactions.
3. Hematoma.

E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician.
2. This device is intended to be used by a trained medical professional.
3. Do not re-sterilize this device.
4. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
5. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.
6. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
7. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
8. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

H. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The contents are manufactured from ABS (acrylonitrile butadiene styrene), Polycarbonate, Polyethylene, Polypropylene, Polyurethane, PVC (polyvinyl chloride), Stainless Steel, and Isoprene Rubber.

J. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

K. INFORMATION

1. **In EU only:** Procedures carried out using these devices may be used on the general population.
2. **In EU only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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