<p><b><i>Symbols Used on Labeling</i></b>  Auf der Verpackung verwendete Symbole  Symboles utilisés sur l'étiquette  Simboli usati sull'etichetta  Símbolos utilizados en las etiquetas  Símbolos utilizados na embalagem</p>
---

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<b>A. INDICATIONS</b> The Arthrex Metal Plate and Screw System is intended to be used for internal bone fixation for bone fractures, fusions, or osteotomies in the ankle, foot, hand, wrist clavicle, scapula, olecranon, humerus, radius, ulna, tibia, patella, calcaneous, femur and fibula. <div style="text-align: right;">Please see identification on label or surgical technique for intended surgical site.</div>
<b>B. CONTRAINDICATIONS</b> 1. Insufficient quantity or quality of bone. 2. Blood supply limitations and previous infections, which may retard healing. 3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation. 4. Any active infection or blood supply limitations. 5. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period. 6. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate. 7. Do not use for surgeries other than those indicated.

<b>C. ADVERSE EFFECTS</b> 1. Infections, both deep and superficial. 2. Foreign body reactions. 3. Injury to nerves, blood vessels, muscles, tendons. 4. Malunion/nonunion of the bone. 5. Necrosis of the region to be treated. 6. Implant failure/breakage. 7. Risk of chronic pain.
--

<b>D. WARNINGS</b> 1. All metallic implant devices used for this surgical procedure should have the same metallurgical composition. 2. Postoperatively and until bone healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device. 3. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device. 4. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Removal of the device should be performed using standard surgical practices for device removal. Implant removal should be followed by adequate postoperative management. 5. The appropriate Arthrex delivery system is required for proper insertion of the device. 6. Detailed instructions on the use and limitations of this device should be given to the patient. 7. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user. An internal fixation device must never be re-used. 8. This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating or migration in the MR environment. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging. 9. Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
--

<b>E. PRECAUTIONS</b> 1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration. 2. Use the appropriately sized drill bit for the screws. 3. Damage to the driver or the screw may result from failure to seat the driver fully into the screw or to align the driver properly with the screw. 4. Do not bend the plate near the locking hole. Bending the plate near the locking hole can distort the holes threading, which prohibits insertion of the screw. 5. Repeated bending of the plate at the same location, or by creating excessive acute angles may potentially lead to premature plate fatigue, failure and or breakage in situ. 6. Prior to drilling the hole depth and diameter has to be identified and appropriate screw needs to be chosen to prevent damages on surrounding tissue. 7. Screws should be inserted by hand and not with powered equipment.
---

<b>F. PACKAGING AND LABELING</b> 1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact. 2. Contact Customer Service if the package has been opened or altered.
---

<b>G. STERILIZATION</b> This device may be provided either sterile or non-sterile. Check the package labeling for more information. <div style="text-align: right;">Certain Arthrex devices that may be used during this procedure are provided nonsterile and must be adequately cleaned and sterilized prior to use or re-use. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions. Cooling – The instrument must be adequately cooled, after being removed from the sterilizer. It should not be touched during the cooling process. Do not place the instrument on a cold surface or immerse in a cold fluid.</div>
--

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:			
Gravity-Displacement Steam Sterilization Cycle	Exposure Temperature	Exposure Time	Drying Time
	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C (270°F)	15 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:			
Gravity-Displacement Steam Sterilization Cycle	Exposure Temperature	Exposure Time	Drying Time
	132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

**G1. PACKAGING**  
Single: A standard packaging material may be used. Ensure that the pack is large enough to contain the device without stressing the seals.  
Sets: Devices may be loaded into dedicated trays or general-purpose sterilization trays. Ensure that sharp edges are protected and do not exceed 8.5 Kg/18 lb per tray. Wrap the trays using the appropriate method.

#### CLEANING AND DISINFECTION OF NON-STERILE IMPLANTS

**MANUAL CLEANING**  
1. Use an enzymatic or alkaline cleaning detergent for rinsing and cleaning. Cleaning solutions can include, but are not limited to: ENZOL® enzymatic, Neodisher® Mediclean Forte, and Thermosept® alk clean. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts.**  
2. Rinse the device thoroughly with water following the cleaning process.

**ULTRASONIC CLEANING**  
1. The devices should be placed in an ultrasonic cleaning unit for a minimum of 20 minutes and processed according to the ultrasonic unit's directions.  
2. The device should be rinsed thoroughly with water following the ultrasonic process.

#### AUTOMATIC WASHING

- Load the devices in the washer such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain (fingers should be open and cannulations/holes positioned to drain).
- Run the automatic wash cycle—minimum cycle parameters:
  - 2 minute cold prewash at 68 ± 9°F (20 ± 5°C),
  - 3 minute cleaning wash (enzymatic or alkaline agent) at 140 ± 9°F (60 ± 5°C),
  - 15 second rinse at 140 ± 9°F (60 ± 5°C),
  - 1 minute thermal rinse at 176 ± 9°F (80 ± 5°C),
  - 6 minute drying phase at high temperature.
- Automatic wash cleaning solutions can include, but are not limited to: ENZOL® enzymatic, Neodisher® Mediclean Forte, and Thermosept® alk clean. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Ultem™, Lexan™, and Cyclocac™.**

#### MANUAL DISINFECTION

- Non-sterile implants should be cleaned before disinfection.
- Immerse non-sterile implant in disinfection solutions for a minimum of 20 minutes.
- Suitable disinfection solutions can include, but are not limited to: CIDEX®, WVICIDE®-01, Gigasept®, Kohrosol®, and equivalent products). Use the supplier's instructions for preparing the solution. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Ultem™, Lexan™, and Cyclocac™.**
- After disinfection, the implant should be rinsed with distilled water or preferably demineralized sterile water.
- Dry implant thoroughly.

#### H. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Stainless Steel (Stainless Steel)  
Titanium Alloy (Ti6Al4V)  
Titanium Niobium Nitride Coating (TiNbN)

#### I. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.  
Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile device is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

#### J. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

<b>Deutsch</b>	<b>DFU-0180r0</b>
----------------	-------------------

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