English

a maximum temperature rise of 3 °C after 15-minutes of continuous scanning.

a. Artifact Information

In non-clinical testing, the image artifact caused by the USP Revision Components extends approximately 60 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. FACTORS AND RISKS IMPACTING THE SAFETY AND SERVICE LIFE OF THE IMPLANT

- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformity of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.
- Thromboses and pulmonary infarction caused during implant site preparation and subsequent / treatment.

H. PACKAGING AND LABELING

- Arthrex implants should be accepted only if the factory packaging and labeling arnt intact.
- All prosthesis parts and components should be stored unopened in the original packaging, and should be kept in protective packaging until immediately prior to use.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbolsglossary**.

I. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

J. MATERIAL SPECIFICATIONS

The revision devices consist of Trunion kits made of titanium alloy, and heads manufactured of cobalt-chromium (Co-Cr) alloy.

K. STORAGE CONDITIONS

Joint prostheses must be treated with great care prior to implantation. Scratches or nicks on the joint surface can result in excessive wear and complications.

- Products must be stored in the original unopened packaging in a dry place and should not be used after the expiration date.
- Do not remove the device from the protective packaging until immediately prior to use.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

M. DIRECTIONS FOR USE

Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.

- The operation should be planned using pre-operative radiographic images. Selection of USP Revision Components must be appropriate to the implants being revised – see guidelines below.
 - If the native glenoid needs to be replaced with an implant or the previous implant needs to be revised, the USP Revision Component head should also be replaced.
 - If the Trunion shows any visible damage, or has become loose after the prosthetic humeral head has been removed, the Trunion should be replaced. If removal of the entire humeral stem is required, replacing with a new Unvers II or Unvers Revers stem is suggested.
- The specific Arthrex implantation instruments are necessary for a revision. The Deflection Torque Driver (3 Nm) must be used for a revision of the Trunion; humeral heads should be extracted with an Arthrex ECLIPSE™ & USP II Extractor. The new humeral head should be fixed with a UNIVERS and ECLIPSE Head Impactor.
- While all heads can be combined with all Trunion kits, the Trunion kits must be selected in accordance with the table below:

Trunion Kit	Long description including the compatible stems	Compatible to the following implanted stems:
AR-9100-20	UNIVERS Trunion Kit 1	USP AR-9100-06 – AR-9100-12 with lot numbers 00.000 - 01.082
AR-9100-21	UNIVERS Trunion Kit 2	USP AR-9100-06 – AR-9100-12 with lot numbers 01.096 - 02.155
AR-9100-22	UNIVERS Trunion Kit 3	USP AR-9100-06 – AR-9100-09 with lot numbers 03.029 - 03.063
AR-9100-23	UNIVERS Trunion Kit 4	USP AR-9100-10 – AR-9100-12 with lot numbers 03.044 - 03.046
AR-9100-24	UNIVERS Trunion Kit 5	USP AR-9100-06 – AR-9100-09 with lot numbers 03.102 - …
AR-9100-25	UNIVERS Trunion Kit 6	USP AR-9100-10 – AR-9100-15 with lot numbers 03.105 - …

F. MRI SAFETY INFORMATION

I. MR CONDITIONAL

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the USP Revision Components are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the USP Revision Components is expected to produce

allow adjustment for offset and version after implantation (loosening of the set screw greater than ¼ turn may cause potential decoupling of the Trunion from the stem). The Torque Driver is used to lock the set screw located on the lateral aspect of the Humeral Stem. Confirm that the 3 Nm mark on the Torque Driver is rotated past the indicator line to ensure that the set screw is properly tightened. Impact the final Head onto the Humeral Stem using the Unvers and ECLIPSE Head Impactor.

NOTES:

- The Trunion of the USP Revision Components must be fully seated on the resected bone.
- The screws which lock the indication block and Trunion must be tightened using the Arthrex Torque Driver, found within the Unvers II instrument set.
- Failure to achieve the appropriate torque requirements when tightening the locking screws may result in the premature loosening of the device.
- Proper locking of components is essential for enduring fixation of the prosthesis.
- The batch number of the implants used must be documented in the corresponding patient records for traceability, together with labels accompanying the implant.
- While unpacking the prosthesis, check that the implant corresponds to the designation on the packaging, and ensure that the implant is treated in a proper, qualified manner.
- The contents of the package should not be implanted after its expiration date.

N. POST-OPERATIVE PROTOCOL

- Patients who have undergone arthroplasty should be routinely examined by the physician to monitor healing, and potential complications.
- The patient must be instructed to immediately inform the surgeon if any adverse changes to the operative joint. Additionally, any falls or other impacts to the operative site should be reported, so that the surgeon can assess whether an additional examination is required. At a minimum, a control X-ray should be sent to the surgeon every 12 months for a radiologic assessment of the operative joint.

- Patients who have undergone arthroplasty should be routinely examined by the physician to monitor healing, and potential complications.
- The patient must be instructed to immediately inform the surgeon if any adverse changes to the operative joint. Additionally, any falls or other impacts to the operative site should be reported, so that the surgeon can assess whether an additional examination is required. At a minimum, a control X-ray should be sent to the surgeon every 12 months for a radiologic assessment of the operative joint.

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