English

A. DEVICE DESCRIPTION

The Bone Marrow Aspirate (BMA) Kit contains a bone marrow aspiration needle, a syringe and a preparation tray.

B. INDICATIONS

The Bone Marrow Aspirate (BMA) Kit is indicated for the aspiration of bone marrow, autologous blood, plasma or other blood components with or without pre-filling with a bone-void filler (allograft, autograft, or synthetic bone graft material).

The Kit is designed to facilitate the pre-mixing of autologous blood, bone marrow aspirate, I.V. fluids, platelet-rich plasma, or plasma with bone-void filler and the delivery to an orthopedic surgical site, as deemed necessary by the clinical use requirements.

C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.
3. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
4. Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Foreign body reactions.
3. Hematoma.

E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician.
2. Procedures carried out using these devices may be used on the general population.
3. The clinical benefits associated with the use of these devices outweigh the known clinical risks.
4. There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
5. This device is intended to be used by a trained medical professional.
6. Do not re-sterilize this device or any of its components.
7. Do not reuse this device or any of its components.
8. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
9. Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
10. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised

mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

11. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
12. Detailed instructions on the use and limitations of this device should be given to the patient.
13. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
14. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. PRECAUTIONS

1. Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

H. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.
Needle: Surgical-grade stainless steel.
Syringe: Polypropylene with a polyisoprene seal.

Preparation Tray: Polypropylene.

J. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

K. DIRECTIONS FOR USE

Users of this device are encouraged to contact their Arthrex representatives if, in their

professional judgment, they require more comprehensive aspiration procedures.

1. Basic Aspiration Procedure

- a. Locate and prepare the bone marrow entry site.
- b. Insert the aspiration needle.
- c. Attach the male connector of the syringe to the female connector handle of the aspiration needle.
- d. Aspirate.
- e. Transfer from the syringe as deemed necessary.

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