English

A. DEVICE DESCRIPTION

The cobalt chrome Univers™ CA (Cuff Arthropathy) head is designed with a larger area of articulation to allow articulation with the acromion in patients with gross rotator cuff deficiency. The CA head is available in many nominal sizes, and is designed to be used with the Univers Shoulder Prosthesis.

B. INDICATIONS

The *Univers CA Heads* are indicated for use in hemi-shoulder replacement. Hemi-shoulder replacement is indicated for:

- A severely painful and/or disabled joint resulting from osteoarthritis, traumatic arthritis or rheumatoid arthritis;
- Fracture-dislocations of the proximal humerus where the articular surface is severely comminuted, separated from its blood supply or where the surgeon’s experience indicates that alternative methods of treatment are unsatisfactory;

Other difficult clinical problems where shoulder arthrodesis or resection arthroplasty are not acceptable (e.g. revision of a failed primary component).

Hemi-shoulder replacement is also indicated for:

- Ununited humeral head fractures
- Avascular necrosis of the humeral head
- Rotator cuff tear arthropathy
- Deformity and/or limited motion

The Univers CA Heads are intended to be used with the Univers Shoulder Stem for cemented or uncemented fixation.

The Univers CA Heads are intended for hemiarthroplasty in shoulder arthritis and trauma where the rotator cuff is not functional and not reconstructable. The CA heads are not intended to be used with glenoid components.

C. CONTRAINDICATIONS

- Patients who are not anatomically and structurally suited to receive the selected implant, or who do not have a functional deltoid muscle.
- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Allergies and other reactions to device materials.
- Loosening of the implant as a result of changed conditions in load transfer, respectively, fatigue wear and breakage of the cement bed and / or tissue reaction to the implant. Loosening is frequently a consequence of one or several of the above listed risk conditions, but can also be caused by inadequate anchoring technique (see #4 below).
- Dislocation, subluxation, or inadequate scope of movement as a result of failure to achieve optimum positioning of the implant.
- Bone fractures as a result of one-sided overload or weakened bone substance.
- Temporary or permanent nerve damage as a result of pressure or hematoma.
- Cardiovascular diseases including venous thrombosis, pulmonary embolism, and cardiac arrest.
- Wound hematoma and delayed wound healing.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be re-used. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Do not re-sterilize this device.
- Detailed instructions on the use and limitations of this device, the patient leaflet and the patient implant card should be given to the patient.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen

prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.

- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Removal of the device should be performed using standard surgical practices for device removal.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The following operative situations may cause premature loosening and complications:
 - Extreme weakening of the bone structure in preparing the bone bed;
 - Unsuitable selection of the implant size;
 - Inadequate cleaning of the bone bed prior to implantation; and,
 - Excessive use of force in placing or fastening the implant, provoking splintering fractures, or causing the bone to tear.
- The appropriate Arthrex delivery system is required for proper insertion of the implant.
- Only Arthrex delivery systems, instruments, and trial prostheses should be used for the implantation procedure.
- Endoprostheses may not be processed mechanically or changed in any other way.
- Do not implant any parts that have been scratched or damaged. An artificial joint is subject to wear and/or can loosen over a period of time. Wear and loosening may make it necessary to re-operate on an artificial joint.
- An infection in an artificial joint may lead to implant removal.
- This device should only be used in conjunction with other implants designed specifically for use with this system.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Univers CA Heads is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Arthrex Univers CA Heads is expected to produce a maximum temperature rise of 3°C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Arthrex Univers CA Heads extends approximately 60 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformation of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.
- Thromboses and pulmonary infarction caused during implant site preparation and subsequent /treatment.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.

- Contact customer service if the package has been opened or altered.
 - All prosthesis parts and components should be stored unopened in the original packaging, and should be kept in protective packaging until immediately prior to use.
 - All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.
- #### I. STERILIZATION
- This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.
- Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.
- #### J. MATERIAL SPECIFICATIONS
- Refer to the package label for the materials.
- The humeral device is made of cobalt-chromium alloy.
- #### K. STORAGE CONDITIONS
- Joint prosthesis must be treated with great care prior to implantation. Scratches or nicks on the joint surface of the prostheses can result in excessive wear and complications.
- Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
 - Do not remove the device from the package until immediately before use.

L. INFORMATION

- In EU Only:** Procedures carried out using these devices may be used on the general population.
- In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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