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All symbols shown may not apply to this device. Please refer to the package labeling for utilized symbols.

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
Arthrex®
Gemini™ SR8

Self-Retaining 8.25mm Cannula System
with No-Squirt Cap

AR-6572






Manufactured by **FRANTZ * MEDICAL**

Patents US 8,360,969, others pending

Please read all information carefully.
Failure to follow the instructions may lead to user and/or patient injury.

NOTE: This insert is intended to provide instructions for assembly and use of the Gemini Cannula System. It is not a reference for surgical techniques.

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Prodotto da:

Frantz Medical Development **CE 0470**
7740 Metric Drive
Mentor, OH 44060, USA
+1 440-255-4514

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Uncaptured wings are UNACCEPTABLE

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ENGLISH

Indications

The Gemini Cannula System has application in arthroscopic orthopedic procedures to establish a path to pass minimally invasive instruments and/or fluid into and out of the surgical site.

Contraindications

This device is not intended for use when minimally invasive techniques are contraindicated.

Device Description

The Gemini Cannula System is a sterile, single patient use instrument consisting of an inner and outer cannula assembly and obturator. The obturator has a blunt plastic tip as to not damage bone, ligament or soft tissue during insertion. The outer cannula has two wings at the distal end that may be deployed to act as retention members, or retracted to allow cannula release, as needed. The inner cannula is capable of telescoping beyond the distal end of the outer cannula to lend the desired position within the surgical site. The cannula pathway permits instrumentation to be passed to the surgical site from outside the patient's body. The cannula pathway contains seals that accommodate instruments up to 8 mm in diameter, and limit fluid leakage when instruments are inserted and withdrawn. A luer port and barb are compatible with standard fittings and provide attachment for fluid inflow and outflow.

DEVICE ASSEMBLY INSTRUCTIONS

Obturator Insertion

Using sterile technique, remove the components from the package.

While grasping the *Inner Cannula only*, verify that the V-notch in the Locking Release Button ① is aligned with the red Inner Cannula triangle ②. Align the red Obturator arrow ③ with the V-notch and triangle, and insert the Obturator into the cannula assembly through the instrument seals.

Wing Capture

Alignment of the Locking Button V-notch, Inner Cannula triangle and Obturator arrow should ensure wing capture. However, before the Obturator is fully inserted through the Instrument Seals, verify that the Outer Cannula wings are upright for capture by the Obturator. If necessary, pinch the wings together ④ so that they are parallel to each other and the center axis of the cannula.

Fully insert the Obturator until the handle of the Obturator contacts the Inner Cannula ⑤ and the wings are fully captured ⑥.

Verification of Proper Assembly

Ensure the wings are FULLY captured by the distal tip of the Obturator.

Secure the Luer Cap onto the Luer Port if the portal is not intended for use as an inflow or outflow pathway.

DEVICE IS READY FOR INSERTION

INSTRUCTIONS FOR USE

Assembly Insertion ⑦

After an entry incision has been made, insert the Inner and Outer Cannula assembly into the joint. If a cannulated Obturator is utilized, the assembly may be placed over a Switching Stick for advancement into the surgical site.

NOTE: The cannulated Obturator is to be used with the 2.6 mm Gemini Switching Stick, or other Switching Stick of comparable or smaller diameter.

Obturator Removal and Wing Deployment ⑧

Once the cannula assembly is in the correct position within the joint, hold the Inner Cannula to retain cannula position and remove the Obturator from the assembly.

The Obturator must be removed before wing deployment. To deploy the wings, use a syringe injection motion and push the Inner Cannula forward with the palm of the hand while pulling against the Outer Cannula towards the palm of the hand with the index and middle fingers.

NOTE: If excessive force is required for wing deployment, verify that the cannula is properly placed within the surgical site and that the Obturator has been removed before continuing to apply force.

For Use as Portal ⑨

The cannula system is now ready to be used as a portal for passing instrumentation and fluid into and out of the joint.

At this time, the No-Squirt Cap supplied with the cannula system may be fitted over the Inner Cannula to mitigate inadvertent spray from the pressurized joint through the seal area. DO NOT place the No-Squirt Cap prior to Obturator insertion.

Assembly Removal ⑩

To remove the cannula assembly from the joint, depress the Locking Release Button while pulling the Inner Cannula away from the patient. As the Inner Cannula is removed, the wings collapse for easy Outer Cannula abstraction.

Manufactured for: Arthrex, Inc.
1370 Creekside Blvd.
Naples, FL 34108, USA
+1 800-934-4404
www.arthrex.com

European Representative: Arthrex GmbH
 Erwin-Hielscher-Strasse 9
81249 München, Germany
Tel: +49 89 909005-0
www.arthrex.de

Manufactured by:  Frantz Medical Development  0470
7740 Metric Drive
Mentor, OH 44060, USA
(440) 255-4514

©2013 Frantz Medical Development. All rights reserved. 201014B 4/15

Warning and Precautions

Arthroscopic procedures should only be conducted by persons with proper knowledge and training techniques. Prior to a minimally invasive, arthroscopic procedure, it is the responsibility of the person performing the procedure to be familiar with all medical literature, methods, and potential complications and hazards that may be associated with such a procedure.

The cannulated obturator is to be used with the Gemini 2.6 mm switching stick or other switching stick of comparable or smaller diameter.

Be sure the obturator and cannula have been properly assembled. Do not attempt to insert the cannula assembly without the obturator.

Refer to the Instructions for Use, and verify that:

The wings are fully captured in the dovetail of the obturator;

The entire assembly is locked together, and that the wings cannot be freed from their captured position by reasonable torsional or tensile forces.

Applying undo force to the cannula assembly and/or obturator during a procedure may damage the product and/or harm the patient. Wing deployment and inner cannula telescoping should not require excessive force. To minimize the risk of injury be sure to:

Confirm locations of key anatomical features;

Direct the distal end of the cannula system away from major anatomical structures;

Verify proper subcutaneous placement of the cannula before attempting to deploy the wings or telescope the inner cannula;

Confirm that the obturator has been removed prior to wing deployment;

Only apply reasonable, necessary pressure to the device.

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The inner cannula allows instruments up to 8 mm to pass through the portal. An attempt to insert any instrument larger than the recommended diameter may result in the inability to pass the instrument, and/or in seal damage and subsequent excess fluid leakage. Special care should be taken during the introduction or removal of any instrument with sharp or angled edges, so as to not inadvertently damage the seal.

Do not place the no-squirt cap prior to obturator insertion.

After removing part or all of the Gemini cannula system, examine the product for damage. Do not reinsert a damaged product. If a product appears damaged after use, make sure no patient and/or user injury has subsequently occurred, and discard the product.

Instruments or devices that come into contact with bodily fluids may require special disposal handling to prevent biological contamination.

This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness, or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, which may result in injury, illness, or death of the patient.

How Supplied

The Gemini Cannula System is supplied sterile for single patient use. Discard after use.

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Representante europeo: Arthrex GmbH