English

A. DEVICE DESCRIPTION

The Double Syringe Autologous Conditioned Plasma (ACP®) System is comprised of an inner syringe inside an outer syringe. The inner syringe was specifically designed to fit within and function in concert with the specifically designed outer syringe. The two syringes are connected via a luer lock connector. A threaded screw cap connected to the luer tip of the outer syringe can be removed prior to use. The Double Syringe Autologous Conditioned Plasma (ACP) System was designed to be used with the existing Rotofix 32 A centrifuge. To prepare the PRP, a sample of the patient's blood can be combined with Anticoagulant Citrate Dextrose A Solution (ACD-A) to prevent coagulation during processing.

B. INDICATIONS

The Arthrex Double Syringe (ACP) Kit is indicated for the safe and rapid preparation of autologous platelet rich plasma (PRP) from a small sample of peripheral blood at the patient's point of care. The PRP is mixed with autograft and/or allograft bone prior to application to a bony defect for improving handling characteristics

C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Allergies and other reactions to device materials.
3. Hematoma
4. Damage to blood vessels and nerve damage resulting in pain or numbness from autologous sampling.
5. Delayed wound healing.

E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician.
2. This device is intended to be used by a trained medical professional.
3. Do not re-sterilize this device.
4. Do not re-use this device.
5. Use only the Rotofix 32 A Benchtop Centrifuge, and follow the manufacturer's instructions. Outcomes using centrifuges from other manufacturers are unknown.
6. Follow the manufacturer's directions that are provided with the package insert for the Anticoagulant Citrate Dextrose Solution A (ACD-A).
7. The safety and effectiveness of this device for bone healing and hemostasis have not been established.
8. The PRP must be used within four (4) hours of blood collection.
9. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in

irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

10. Detailed instructions on the use and limitations of this device should be given to the patient.
11. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
12. ACD-A anticoagulant and PRP prepared are not for intravenous use. For additional information on ACD-A anticoagulant, refer to the component label.
13. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
14. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. PRECAUTIONS

1. Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

H. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.
Outer Syringe: Polypropylene, TPV, Silicone Oil

Inner Syringe: Polypropylene, Silicone Oil

Cap: Polyethylene
ACD-A: See component label

Note: The ACD-A included in the kit is only for use with Arthrex PRP Systems. Discard the unused portion.

J. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

K. INFORMATION

1. **In EU Only:** Procedures carried out using these devices may be used on the general population.
2. **In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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