## English

### A. DEVICE DESCRIPTION

The Arthrex iBalance® High Tibial Opening Wedge Osteotomy (HTO) System includes instrumentation, implant, and fixation anchors, and a surgical instrument system. The wedge-like implants, and cortical and cancellous fixation anchors are available in various sizes.

### B. INDICATIONS

The Arthrex iBalance High Tibial Opening Wedge Osteotomy (HTO) System is comprised of an anatomically designed implant, fixation anchors, and a surgical instrument system. The wedge-like implants, and cortical and cancellous fixation anchors are available in various sizes.

### C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions which tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.
- Patients who have not reached skeletal maturity.

### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

### E. WARNINGS

- An internal fixation device must never be reused. The devices may be contaminated and not properly sterilized. This may result in infection, toxic reaction or similar adverse effect. Re-use of a device may result in mechanical damage such as deformation of the threads on the implant or anchors, weakening of the anchor, or deformation of the feature that inhibits back-out of the anchor. This may cause difficulties in implantation or failure of the device post-operatively.
- Do not re-sterilize this device.
- Post-operatively, the patient is placed in a post-op, hinged knee brace locked in 0 degrees of extension.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex iBalance instrument system is required for proper implantation of the device. Weight bearing should not begin until there is evidence of bone growth.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2)

Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or break-age, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

### F. MRI SAFETY INFORMATION

#### 1. MR Conditional

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the iBalance High Tibial Opening Wedge Osteotomy System is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorpion rate (SAR) of 1-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the *iBalance High Tibial Opening Wedge Osteotomy System* is expected to produce a maximum temperature rise of up to 6 °C after 15-minutes of continuous scanning.

#### a. Artifact Information

*In non-clinical testing, the image artifact caused by the iBalance High Tibial Opening Wedge Osteotomy System can extend up to approximately 120 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

### G. PRECAUTIONS

- Before opening the implant sterile barrier pouch, confirm implant size by viewing the engraving on the implant through the clear packaging.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

**CAUTION:** Improper alignment of the instrumentation may result in unintended tibial A/P slope change.

### H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbolsglossary**.

### I. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

Certain Arthrex devices and instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities”

### J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The devices are made out of polyetheretherketone (PEEK) and tantalum.

### K. STORAGE CONDITIONS

This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

### L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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## iBalance® High Tibial Opening Wedge Osteotomy System

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
This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the “Arthrex U.S. Product Warranty” section of the Arthrex, Inc. website, found at [www.arthrex.com](http://www.arthrex.com) whose provisions are incorporated herein by reference.

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