# English

#### A. DEVICE DESCRIPTION

The SpeedBridge™ Implant Systems are available in multiple configurations for use in the SpeedBridge technique for the shoulder. The configurations consist of SwiveLock® suture anchors (Bio, BioComposite, PEEK, or Self-Punching) with suture and manual surgical instruments.

The Achilles SpeedBridge Implant Systems are available in multiple configurations for use in the SpeedBridge technique for Achilles repair. The configurations consist of the BioComposite SwiveLock suture anchorswith suture and without sutureand manual surgical instruments.

The Achilles SutureBridge™ Implant System consists of Bio-Corkscrew® and Bio-PushLock™ suture anchors with suture and manual surgical instruments for use in the SutureBridge technique for Achilles repair.

The Eclipse™ Subscapularis SpeedBridge Kit consists of SwiveLock and FiberTak® DR (self-punching) suture anchors with suture and manual surgical instruments.

#### B. INDICATIONS

Arthrex Corkscrew®, PushLock® and SwiveLock suture anchors are intended to be used for soft (soft tissue) fixation to bone in the shoulder, foot and ankle.

FiberTak DR suture anchor is intended for fixation of suture (soft tissue) to bone in the shoulder (rotator cuff repair).

The Arthrex Suture Family is intended for use in soft tissue approximation and or ligation. The sutures may be incorporated, as components, into surgeries where constructs including those with allograft or autograft tissues are used for repair.

#### C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any acute infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.
- Arthrex suture is not for use in cardiac indications.

#### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Adverse reactions have not been noted with Arthrex suture products in animal testing.
- Common non-absorbable suture reactions may include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site.
- Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.
- Silicone sensitivity, though very rare, has been reported.

#### E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be reused.
- Do not re-sterilize this device.
- Attempting implantation into hard, dense bone and/or drilling/punching smaller diameter holes than recommended may cause failure (breakage) of the implant during insertion.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device, and the patient leaflet ([www.arthrex.com/patientleaflets](http://www.arthrex.com/patientleaflets)) and the patient implant card should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

- Once open, discard unused suture.
- Do not expose suture to heat.
- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Arthrex FiberWire® or Arthrex UHMWPE suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.
- As with any foreign body, prolonged contact of suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formations. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.
- The safety and effectiveness of these devices for use as an artificial ligament or tendon has not been established.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

#### F. MRI SAFETY INFORMATION

##### 1. MR Conditional: Self-Punching SwiveLock

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the metal (titanium) portion of the SwiveLock suture anchor is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the SwiveLock suture anchors are expected to produce a maximum temperature rise of up to 1.8 °C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the SwiveLock suture anchor as extend up to approximately 17 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

##### 2. MR Safe

*The Corkscrew, PushLock, and SwiveLock suture anchor devices manufactured from only polyetheretherketone (PEEK) poly (L-Lactide, PLLA), Poly (L-Lactide-co-D acid, PLDLA), and/or Poly (L-Lactide acid, PLLA) and tri-calcium phospate (TCP) are MR safe.*

*The FiberTak DR suture anchor and Arthrex Suture manufactured from UHMWPE and polyester with or without silicone elastomer coating, cyanoacrylate, and nylon are MR safe.*

#### G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate suture anchor size based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Make sure to use the recommended drill bit or punch to create the bone socket.
- PushLock and SwiveLock suture anchors only:** During anchor insertion, ensure that the angle of anchor insertion is coaxial to that of the previously prepared bone socket.
- PushLock and SwiveLock suture anchors only:** Insert the driver into the bone socket until the anchor body makes contact with the bone. Preview and adjust suture tension, if necessary. Tension will not increase during final advancement of the anchor body.
- PushLock and SwiveLock suture anchors only:** Ensure that the anchor body is in full contact with the bone before advancing the anchor body into the prepared bone socket.
- Self-punching suture anchors only:** Insertion in very hard bone may require pre-drilling or pre-punching a bone socket to avoid damage to the implant.
- Self-punching suture anchors only:** Ensure that the angle of anchor insertion is perpendicular to the bone.
- In handling suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- Assure that all knots have been secured using accepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Care should be taken to prevent damage to surrounding tissue or user puncture due to improper handling of the needlepoint.
- Do not grasp any suture needle if present at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

#### H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

#### I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

#### J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

These suture anchor devices consist of either one or two components. Each component is manufactured from either/or titanium alloy, polyetheretherketone (PEEK), Poly (L-Lactide acid, PLLA), Poly (L-Lactide-co-D acid, PLDLA), or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP).

**Suture:** See package label for size and type of suture provided with device.

The FiberWire, FiberTap®, and TigerTape® sutures are made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. The FiberWire and TigerWire® sutures include FiberStick™, TigerStick®, FiberChain®, FiberLink™, TigerLink™, TigerTail®, FiberLoop® and TigerLoop™ sutures. The FiberLink, TigerLink, FiberLoop and TigerLoop sutures may also be braided from SutureTape sutures. Additional materials may include silicone elastomer coating (except suture with the suffix -Tape), cyanoacrylate, and may include nylon.

The sutures supplied meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The suture dyes may include: D&C Blue No. 6, D&C Green No. 6, Logwood Black and FD&C Blue No. 2. Suture strands that are dyed with Logwood black are made of nylon. Other non-absorbable suture is made from polyester, polytetrafluoroethylene (PTFE) coating, and cyanoacrylate.

Arthrex sutures, when tested per ISO 10993, Biological Evaluation of Medical Devices-Part 10: Tests for Irritations and Sensitization, had no reactions of allergic or sensitive nature. The dyed suture and coating are pharmacologically inactive.

Arthrex FiberWire suture is not absorbed, but may become encapsulated in the surrounding connective tissues. Arthrex FiberWire is not known to have significant change in tensile strength in vivo.

#### K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C /90° F, and should not be used after the expiration date.

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

#### L. INFORMATION

- In EU Only:** Procedures carried out using these devices may be used on the general population.
- In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

Arthrex FiberWire suture is not absorbed, but may become encapsulated in the surrounding connective tissues. Arthrex FiberWire is not known to have significant change in tensile strength in vivo.

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

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