# English

### A. DEVICE DESCRIPTION

The *InternalBrace*™ Ligament Augmentation Implant Systems are available in multiple configurations for use in ligament repair. The *InternalBrace* Ligament Augmentation Repair Implant System consists of a BioComposite SwiveLock® anchor, a BioComposite SwiveLock anchor with TigerTape™, and manual surgical instruments. The Hand/Wrist *InternalBrace* Implant System consists of PEEK SwiveLock SL anchors, suture, and manual surgical instruents. The Knee Ligament *InternalBrace* Implant System consists of BioComposite SwiveLock anchors, suture, and manual surgical instruments. The Forefoot *InternalBrace* Implant Systems consist of Tenodesis screws (PEEK or BioComposite), a button, suture, and manual surgical instruments.

### B. INDICATIONS

The SwiveLock suture anchors are intended to be used for suture (soft tissue) fixation to bone in the foot, ankle, knee, hand, and wrist. The Arthrex Suture Family is intended for use in soft tissue approximation and/or ligation. The sutures may be incorporated, as components, into surgeries where constructs including those with allograft or autograf tissues are used for repair.

The Tenodesis Screw is intended to provide soft tissue reattachment (i.e. fixation of ligament and tendon graft tissue).

The Suture Button is used for fixation of bone to bone or soft tissue to bone, and is intended as fixation posts, a distribution bridge, or for distributing suture tension over areas of ligament or tendon repair.

### C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.
- Arthrex suture is not for use in cardiac indications.

### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Adverse reactions have not been noted with the Arthrex FiberWire® product in animal testing.
- Common non-absorbable suture reactions may include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site.
- Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.
- Sutures Only:** Silicone sensitivity, though very rare, has been reported.

### E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- Attempting implantation into hard, dense bone and/or drilling/punching smaller diameter holes than recommended may cause failure (breakage) of the implant during insertion.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device, the patient leaflet and the patient implant card should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
- Once open, discard unused suture.

- Do not expose suture to heat.
- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Arthrex FiberWire or Arthrex UHMWPE suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.
- As with any foreign body, prolonged contact of suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formations. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.
- The safety and effectiveness of these devices for use as an artificial ligament or tendon has not been established.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

### F. MRI SAFETY INFORMATION

#### 1. MR Conditional

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Suture Buttons manufactured from Stainless Steel are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Suture Buttons are expected to produce a maximum temperature rise of up to ≤2.5 °C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the Suture Buttons can extend up to approximately 17 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

#### 2. MR Safe

*The SwiveLock suture ancbor devices manufactured from only polyetheretherketone (PEEK) poly (L-Lactide, PLLA), Poly (L-Lactide-co-D acid, PLDLA), and/or Poly (L-Lactide acid, PLLA) and tri-calcium phsphate (TCP) are MR safe.*

*Arthrex Suture manufactured from UHMWPE and polyester with or without silicone elastomer coating, cyanoacrylate, and nylon are MR safe.*

*The Tenodesis screws manufactured from only poly-etheretherketone (PEEK), polypropylene or poly (L-Lactide), poly (L, DL-lactide; PLDLA), biphasic calcium phsphate, biphasic tricalcium phsphate (TCP) and/or hydroxyapatite (HA) are MR safe.*

### G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate suture anchor size based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Make sure to use the recommended drill bit or punch to create the bone socket.
- During anchor insertion, ensure that the angle of anchor insertion is coaxial to that of the previously prepared bone socket.
- Insert the driver into the bone socket until the anchor body makes contact with the bone. Preview and adjust suture tension, if necessary. Tension will not increase during final advancement of the anchor body.
- Ensure that the anchor body is in full contact with the bone before advancing the anchor body into the prepared bone socket.
- Insertion in very hard bone may require pre-punching a bone socket to avoid damage to the implant.
- Ensure that the angle of anchor insertion is perpendicular to the bone.
- Under insertion of the device may leave the proximal end of the implant protruding beyond the cortical bone, which could potentially cause soft tissue irritation and/or pain post-operatively.
- Bio-Tenodesis™ Screw only:** Use the appropriate size Arthrex Drill to create a pilot hole in the bone.
- In handling suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- Assure that all knots have been secured using accepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting modifications. Care should be taken to prevent damage to surrounding tissue or user puncture due to improper handling of the needlepoint.
- Do not grasp any suture needle if present at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

### H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbology](http://www.arthrex.com/symbology).

### I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

### J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The SwiveLock suture anchor devices consist of either one or two components. Each component is manufactured from either/ or Poly (L-Lactide acid, PLLA), Poly (L-Lactide-co-D acid, PLDLA), or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP).

**Suture:** See package label for size and type of suture provided with device.

FiberWire, FiberTape, and TigerTape sutures are made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. Additional materials may include silicone elastomer coating (except suture with the suffix -Tape), cyanoacrylate, and may include nylon.

The sutures supplied meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The suture dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon. Other non-absorbable suture is made from polyester, polytetrafluoroethylene (PTFE) coating, and cyanoacrylate.

Arthrex sutures, when tested per ISO 10993, Biological Evaluation of Medical Devices-Part 10: Tests for Irritations and Sensitization, had no reactions of allergic or sensitive nature. The dyed suture and coating are pharmacologically inactive.

FiberWire suture is not absorbed, but may become encapsulated in the surrounding connective tissues. FiberWire suture is not known to have significant change in tensile strength in vivo.

**Tenodesis screws** are manufactured from polyetheretherketone (PEEK), polypropylene or poly (L-Lactide). Bio-composite devices are made of made of various configurations of polymers including poly (L-lactide), poly (L, DL-lactide; PLDLA), biphasic calcium phosphate, biphasic tricalcium phosphate (TCP) and hydroxyapatite (HA).

**Button:** This device is made of stainless steel.

### K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-Bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

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