A. DEVICE DESCRIPTION

The Thrombinator™ System for use with Arthrex PRP Systems (Arthrex Angel® concentrated Platelet Rich Plasma (cPRP) System or Arthrex Double Syringe (ACP) System) includes a sterile, single use device and a filter that aids in the preparation of autologous thrombin serum. Additional equipment (i.e., winged infusion withdraw needle, syringes, syringe cap(s), anticoagulant, and wet ice) used for blood drawing and processing is available and sold separately and to be supplied by the end-user in order to operate the device as intended.

B. INDICATIONS

The Thrombinator System for use with Arthrex PRP systems (Arthrex Angel concentrated Platelet Rich Plasma (cPRP) System or Arthrex Double Syringe (ACP) System) is designed for the preparation of autologous serum from anticoagulated or non-anticoagulated peripheral blood, platelet poor plasma, or platelet rich plasma (PRP) that is to be mixed with PRP and autograft or allograft bone prior to application to a bony defect for improving handling characteristics.

C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.
3. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
4. Do not use for surgeries other than those indicated.
5. Use as a dialyzer or for dialysis with a dialysate.
6. Direct connection to patient's vascular system of circulating blood volume.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Allergies and other reactions to device materials.
3. Hematoma.
4. Damage to blood vessels and nerve damage resulting in pain or numbness from autologous sampling.
5. Delayed wound healing.

E. WARNINGS

1. The autologous thrombin serum exhibits thrombin activity and therefore must not be injected into or otherwise allowed to enter the patient's vascular system.
2. Do not re-sterilize this device.
3. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.

4. The safety and effectiveness of the device for bone healing and hemostasis has not been established.

F. MRI SAFETY INFORMATION

This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown.

G. PRECAUTIONS

Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

The device should be operated according to the surgical technique provided with the device.

H. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbols-glossary.

I. STERILIZATION

The device is supplied sterile. Refer to the package label for the sterilization method.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Thrombinator Construct: Polycarbonate, silicone

Filter: Polyethylene, stainless steel, STYRENE-BUTADIENE COPOLYMER, AND POLYCARBONATE

Vent: Acrylic copolymer and nylon

Beads: Borosilicate Glass

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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