English

A. DEVICE DESCRIPTION

The Compression and Small Staples are metal devices.

The Plaple® is a metal plate/screw combination device with one end for metal cortical fixation.

The DynaNite™ Nitinol Staple provides a means of bone fixation for the management of fractures and reconstructive surgery.

- DynaNite Staples are designed to exhibit super-elastic properties.
- The staples are sold in two configurations:
 - 1. as part of a single use disposable kit - the kit includes a staple, pre-loaded onto a delivery device, a drill bit, a drill guide, alignment pins and a tamp and

- as a staple loaded on a delivery device – to be procured labels for the sterilization method.

DynaNite Staples must not be re-sterilized.

These devices are intended to be used for fixation such as: Lisfranc arthrodesis, mono- or bi-cortical osteotomies in the foot/ankle, first metatarsophal arthrodesis, skin osteotomy, midfoot and hindfoot arthrodesis or osteotomies, fixation of osteotomies for hallux valgus treatment (Scarf and Chevron) and arthrodesis of the metatarsocuneiform joint to osteotomy and stabilize metatarsus primus varus.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities for long durations during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.

- Do not use for surgeries other than those indicated.
- DynaNite Staple only:** Comminuted bone surface that would not allow for staple placement.
- DynaNite Staple only:** Pathologic bone conditions such as osteopenia that would impair the ability to securely fix the implant.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
 - Foreign body reactions.
- E. WARNINGS**
- An internal fixation device must never be reused.
 - All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
 - The Arthrex Plaple is intended for use only with 2.3 mm diameter LPS cortical stainless steel screws.
 - Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.

- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.

- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.

The quality of water utilized for the cleaning and disinfection processes should be validated at the end-user site to ensure that residuals will not interfere with subsequent processing steps in accordance with ANSI TR30-2011. Deionized water was the water quality used for manual and automated cleaning process validations performed to establish the required processing parameters as outlined in the "Cleaning" section of this document.

- DynaNite Staple only:** The implants are not designed to replace normal healthy bone or withstand the stress placed upon the device by full or partial weight-bearing or load bearing in the presence of nonunion, delayed union or incomplete healing. Immobilization of the treatment site using routine methods (casting, splinting, etc.) should be maintained until no bone healing has occurred (4-6 weeks).
- DynaNite Staple only:** Reduction of the fracture should be achieved and maintained prior to implanting the device. The staple should not be relied upon to achieve closure or reduction of the fracture line.

- DynaNite Staple only:** Additional processing or reprocessing of the implant, may affect the material properties of the Nitinol, changing or otherwise reducing the effectiveness of the implant.

- DynaNite Staple only:** Drill bits are designed for use specifically with the Arthrex supplied drill guides. Reprocessing may affect the compatibility with other instruments and usability of the instruments.
- DynaNite Staple only:** If sterilization is compromised prior to insertion, a different sterile implant will need to be used. Do not re-sterilize or re-use.

F. MRI SAFETY INFORMATION

I. MR CONDITIONAL.

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Staples are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, only
- Maximum static gradient magnetic field of 2000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

- Under the scan conditions defined, the Arthrex Staples is expected to produce a maximum temperature rise of up to 14.5° C after 15-minutes of continuous scanning.

a. Artifact Information

In non-clinical testing, the image artifact caused by the Arthrex Staples can extend up to approximately 120 cm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides

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II. MACHINE (AUTOMATIC) CLEANING, AND THERMAL DISINFECTION

The instructions for Machine (Automatic) Cleaning and Thermal Disinfection may be followed as an alternative cleaning method to Manual Cleaning and Thermal Disinfection.

- Load the devices in the washer so that all design features of the device are accessible to cleaning and so that design features that might require wash cycle. The minimum cycle parameters are listed below:

Phase	Recirculation Time	Temperature	Detergent
Cold Pre-Wash	3 Minuten	68 ± 9°F (20 ± 5 °C)	N/A
Cleaning Wash	10 Minuten	150 ± 9°F (65.5 ± 5°C)	Enzymatisches oder alkalische agent
Rinse 1	3 Minuten	122 ± 9°F (50 ± 5°C)	N/A
Rinse 2	3 Minuten	122 ± 9°F (50 ± 5°C)	N/A
Thermal Disinfection Rinse	5 Minuten	194 ± 9°F (90 ± 5°C)	N/A
Drying	6 Minuten	239°F (115°C)	N/A

- Automatic wash cleaning solutions can include, but are not limited to: neidisher® Mediclen forte, Thermostop® alkaliclean, Prolysics® Ultra Concentrate Enzymatic Cleaner, and Proklex™ Npl Neutral Detergent. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylenepropylene), ABS (Acrylnitrilbutadiene Styrene), Utem™, Lexan™, and Cytolac™.** If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.
- Check the devices for visible soil. Repeat cleaning if soil is visible and re-inspect.

L. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. This device is made of stainless steel, titanium, Cobalt Chromium or Nitinol.

M. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date. Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited, the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

N. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic forms. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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