English

A. DEVICE DESCRIPTION

The PEEKPower™ Plates are available in different shapes, sizes and orientations (e.g. left and right types). The plates have specific sized holes for screws to provide fixation.

B. INDICATIONS

The PEEKPower Plates are used as an internal fixation device to stabilize fractures or osteotomies at the time of healing. The appropriate plate should be chosen for the site to be treated.

The PEEKPower Humeral Fracture Plate is intended for fixation of proximal fractures of the humerus.

The PEEKPower Osteotomy Plates provide fixation following opening wedge osteotomies.

The plates are intended to be used with adequate post-operative immobilization.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.
- PEEKPower High Tibial Osteotomy (HTO) and Lateral Distal Femoral Osteotomy (LDFO) Plates:** Previous closing wedge osteotomy in the same region of the same extremity.
- PEEKPower HTO and LDFO Plates:** Closing Wedge Osteotomy in the same region of the same extremity.
- PEEKPower HTO and LDFO Plates:** Rotational Osteotomy in the same region of the same extremity.

D. ADVERSE EFFECTS

- Foreign body reactions
- Infections, both deep and superficial.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- Use only the Arthrex Screws provided with the PEEKPower Plate systems.
- Osteotomy plates only:** Corresponding preoperative planning must be followed.
- Postoperatively and until bone healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Removal of the device should be performed using standard surgical practices for device removal. Implant removal should be followed by adequate postoperative management.

- Detailed instructions on the use and limitations of this device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Do not re-sterilize this device.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- The PEEKPower Plates should not be used in patients whose body weight is more than 100 kg (220 lbs.).
- An internal fixation device must never be reused.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the PEEKPower Plates are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the PEEKPower Plates are expected to produce a maximum temperature rise of 5.2 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the PEEKPower Plates extends approximately 66 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- The plate cannot be bent or intra-operatively contoured. Contouring the plate could result in premature failure/breakage.
- The locking screw fixation will be compromised if the plate is damaged during the drilling process.
- The locking screws can be inserted into the same plate hole a second time, at the same or different angle, if necessary. For example, if a locking screw needs to be angled differently in the hole, or a different length screw needs to be used.
- The screws should be inserted by hand with the appropriate screwdriver and not with powered instruments. The most effective fixation is achieved when the screw head is inserted flush with the top surface of the plate.
- Proper placement of the plate is critical to avoid impingement and to ensure optimal screw placement.
- PEEKPower, HTO and LDFO plate: For a fracture of the far cortex, a compression screw must be used to reduce the fracture fragments. Once the fragments are stabilized and reduced, the compression screws must be replaced with locking screws for the most stable fixation.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Surgeons must apply their professional judgment when determining the appropriate plate size based on the specific indication, preferred surgical technique, and patient history.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI S779 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

PEEKPower Plates are made of Carbon Fiber Reinforced PEEK (CF-PEEK) and tantalum.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

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PEEKPower™ Plates

PEEKPower-Platten
Placas PEEKPower
Plaques PEEKPower
Placche PEEKPower

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Tel: +49 89 909005-0

www.arthrex.de

This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the "Arthrex U.S. Product Warranty" section of the Arthrex, Inc. website, found at www.arthrex.com whose provisions are incorporated herein by reference.

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Arthrex GmbH