<p><b>Symbols Used on Labeling</b>  Auf der Verpackung verwendete Symbole  Symboles utilisés sur l'étiquette  Simboli usati sull'etichetta  Símbolos utilizados en las etiquetas  Símbolos utilizados na embalagem</p>
--

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## English DFU-0188r3

**A. DEVICE DESCRIPTION**  
The 90° Hook Electrode is to be used in arthroscopic surgery with an Electrosurgical Pencil and an Arthrex generator.

This device is to be used as a hand-switching monopolar accessory in conjunction with a compatible electrosurgical generator for tissue cutting and coagulation. The Arthrex Hook Electrode has a standard diameter shaft and fits Arthrex, and most other electrosurgical pencils.

**B. INDICATIONS**  
The Arthrex electrosurgical cutting and coagulation devices and accessories are monopolar in structure, allowing for use in a range of surgical procedures. The Arthrex electrosurgical ablation devices and accessories are intended for use in resection, ablation, and coagulation of soft tissue and hemostasis of blood vessels and tissue in general, arthroscopic, and orthopedic procedures.  
Specifically, these devices and their accessories will be used for general surgeries, and open and arthroscopic surgery of the shoulder, wrist, hand, elbow, hip, knee, and ankle.

**C. CONTRAINDICATIONS**  
1. Insufficient quantity or quality of bone.  
2. Blood supply limitations and previous infections, which may retard healing.  
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.  
4. Any acute infection or blood supply limitations.  
5. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.  
6. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.  
7. Do not use for surgeries other than those indicated.

**D. WARNINGS**  
1. This device has been specifically designed for use in arthroscopy while in the presence of conductive fluid (i.e. saline, Lactated Ringer's, etc). Do not use for other procedures.  
2. **Danger - Explosion Hazard:** Do not use in the presence of flammable anesthetics.  
3. **Fire Hazard:** Both oxygen (O<sub>2</sub>) and nitrous oxide (N<sub>2</sub>O) support combustion. Avoid O<sub>2</sub> and N<sub>2</sub>O enriched atmospheres. Enriched atmospheres may result in fires and burns to patients or surgical personnel.  
4. The metal shaft of the electrode must fit completely and securely into the pencil. If metal is visible, do not use the electrode. Exposed metal may result in electric shock or burn to the patient or surgical personnel.  
5. Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings to achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power setting adjustments.

6. **Fire/Explosion Hazard:** The following substances contribute to increased fire and explosion hazards in the operating room:  
• Oxygen enriched environments.  
• Oxidizing agents such as nitrous oxide (N<sub>2</sub>O) atmospheres.  
◦ Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.  
◦ Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.  
◦ If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube.  
◦ If possible, stop supplemental oxygen at least one minute before and during uses of electrosurgery.  
• Alcohol-based skin prepping agents and tinctures.  
• Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.  
• Naturally occurring flammable gases (such as methane) that may accumulate in body cavities.

7. **Fire Hazard:** Always place the active electrode in a clean, dry, insulated safety holster when not in use.  
• Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.  
• Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes. Place longer electrodes such as extended electrodes away from the patient and drapes.

8. **Fire Hazard:** The sparking and heating associated with electrosurgery can provide an ignition source.  
9. Observe fire precautions at all times:  
• When using electrosurgery in the presence of gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes. Tent drapes to allow vapors to mix with room air before using electrosurgery.  
• Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen enriched environments such as in throat or mouth procedures. With sufficient heating, eschar can become a glowing ember and pose a fire hazard both as an ignition source and as a fuel. Keep the electrode clean and free of all debris.  
10. Facial and other body hair is flammable.  
• Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

11. Conductive fluids (e.g., blood or saline) in direct contact with an active electrode or in close proximity to any active accessory may carry electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling with the active electrode or capacitive coupling between the active electrode and the external surface of the electrode insulation. Therefore, to prevent unintended burns in the presence of conductive fluids:  
• Always keep the external surface of the active electrode away from adjacent tissue while activating the electrosurgical generator.  
12. Buzzing the Hemostat or other instrument may cause severe burns to the surgeon or patient, including burns to the surgeon's hands. This practice is not recommended. Some surgeons may disregard this warning and elect to buzz the hemostat during surgical procedures. The hazards of such a practice cannot be eliminated. The following precautions may slightly reduce the risks of burn:  
• Do not "buzz the hemostat" with a Needle electrode.  
• Do not lean on the patient, the table, or the retractors while buzzing the hemostat.  
• Activate cut rather than coag. Cut has a lower voltage that coag.  
• Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.  
• Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.  
• Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes current concentration at the finger tips.  
• "Buzz the hemostat" below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon's hands.  
• When using a coated or nonstick blade electrode, place the edge of the electrode against the hemostat or other metal instrument.

13. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.  
14. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are

not limited to cross-infection, breakage resulting in irretreivable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

15. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.

**E. PRECAUTIONS**  
1. Arthrex recommends maintaining a fluid outflow temperature of no greater than 40°C (104°F) from the surgical site while ablating target tissue.  
2. Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.  
3. Do not activate electrodes while in contact with or near other instruments, including cannulae. Localized burns to the patient or physician may result.  
4. Do not modify or add to the insulation of active electrodes.  
5. Activate the electrosurgical unit only when you are ready to deliver electrosurgical current and the active tip is in view and near target tissue (especially if looking through an endoscope).  
6. Deactivate the electrosurgical unit before the tip leaves the surgical site.  
7. **DO NOT RESTERILIZE.**  
**IMPORTANT:** Clean the electrode often with moist gauze or other material appropriate for the electrode type.

**F. PACKAGING AND LABELING**  
1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.  
2. Contact Customer Service if the package has been opened or altered.

**G. STERILIZATION**  
This device is provided sterile. Refer to the package label for the sterilization method.  
Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

**H. STORAGE CONDITIONS**  
Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

**I. INFORMATION**  
Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

### J. DIRECTIONS FOR USE

**Install the Electrode**  
1. Ensure the pencil is not connected to the generator, or the generator is Off or in the Standby mode, if available.  
2. Grasp the insulating sleeve on the electrode. Insert the electrode into the pencil.  
3. Ensure the electrode is fully inserted into the pencil. The shank and insulating sleeve should fit securely into the active accessory. If the shank and/or insulating sleeve does not fit, or the insulation will not insert 3.2 mm (1/8 in.), do not use this electrode/accessory combination.  
4. Refer to the generator instruction manual for proper procedures for connecting the active accessory to the generator.  
5. A tip protector covers the end of some electrodes. If a tip protector is present, remove and discard it before use. Do not use it for insulation purposes.

**Maximum Power**  
**NOTICE:**  
Exceeding maximum power settings may result in patient injury or product damage.  
Arthrex recommends the following maximum power settings:

Mode	Arthrex Generator
Coag	15 W
Cut	20 W

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