## English

### A. DEVICE DESCRIPTION

Plates are available in different shapes, sizes and orientations (e.g. left and right types). The plates have specified sizes (holes for screws to be used with ANS/AMI 179) with multiple holes (i.e. Arthrex Mesh Plates) may be cut to size.

The Distal Clavicle Plate Button is designed to fit securely in the holes of the fracture plates. The Plate Button is available in one size only. The Distal Clavicle Plate Button is available as a single button or assembled in a construct with Arthrex suture. Washers are available in diameters of 6 – 13 mm.

### B. INDICATIONS

The *Arthrex Fracture Plates and Ankle Plates* Fracture Plates are used for fixation of fractures of the radius, ulna, forearm, wrists, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus and fibula. The Arthrex Fracture Plates are to be used with the *Arthrex Low Profile Screws*.

The *Distal Clavicle Plate Button* is intended for use with the clavicle plates for clavicle indications such as for the treatment of syndesmotic trauma, such as fixation of acromioclavicular separation due to coracoclavicular ligament disruption, and this button may not be used alone. The button is intended to be used with #5 FiberWire® or FiberTape®.

The *Arthrex Low Profile Plates, Arthrex Compression Plates and Calcaneal Fracture Plates* are intended to be used for internal fixation of fractures, fissures, osteotomies and non-unions in the ankle, foot, hand, and wrist, such as opening wedge osteotomies of Hallux Valgus.

The *Arthrex Distal Extremity Plate System* and *Arthrex Mesh Plates* are intended for use in stabilization of fresh fractures, revision procedures, osteotomies, joint fusion and reconstruction of small bones and bone fragments of the hand/ wrist, foot/ankle and osteotopic bone.

The *Arthrex Distal Radius System* is intended for fixation of fractures, fissures, osteotomies and non-unions of the radius, ulna, olecranon, metacarpal, metatarsal and malleolus.

The *Arthrex* washers and certain plates are intended as adjunct devices.

### C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, distort or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

### E. WARNINGS

- An internal fixation device must never be reused.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be taken to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

### F. MRI SAFETY INFORMATION

#### I. MR CONDITIONAL

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Plates are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system.
- Under the scan conditions defined, the Arthrex Plates are expected to produce a maximum temperature rise up to 14.9° C after 15-minutes of continuous scanning.

#### a. Artifact Information

*In non-clinical testing, the image artifact caused by the Arthrex Plates can extend up to 120 mm from the implant. From this point when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

### G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
  - Do not bend the plate near the locking hole. Bending the plate near the locking hole can distort the holes threading, which prohibits insertion of the screw.
  - Repeated bending of the plate at the same location, or by creating excessive acute angles may potentially lead to premature plate fatigue, failure and or breakage in situ.
  - Screws should be inserted by hand and not with powered equipment.
- H. PACKAGING AND LABELING**
- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
  - Contact Customer Service if the package has been opened or altered.
  - All of the symbols used on the labeling along with the title, device designation number may be found on our website at [www.arthrex.com/symbologyssary](http://www.arthrex.com/symbologyssary).

### I. VALIDATION

The instrument cleaning, disinfection, and sterilization methods in this DFU have been validated in compliance with federal and

This device (i.e. for open but unused) may be re-sterilized following cleaning, disinfection prior to use. Sterilization wrap that has been FDA-cleared for the specified sterilization cycles should be used. Ensure that the wrap is large enough to contain the device without stressing the seals.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to the Arthrex Instruments Directions for Use and ANS/AMI ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities” for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
Gravily-Displacement Steam Sterilization Cycle	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C (270°F)	15 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	135°C (275°F)	10 Minutes	30 Minutes
	132°C (270°F)	4 Minutes	20 to 30 Minutes
Pre-vacuum Cycle	135°C (275°F)	3 Minutes	16 Minutes

### STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:

	Exposure Temperature	Exposure Time	Drying Time
Gravily-Displacement Steam Sterilization Cycle	132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C – 135°C (270°F – 275°F)	4 Minutes	20 to 30 Minutes

Cooling – The device must be adequately cooled, after being removed from the sterilizer. It should not be touched during the cooling process. Do not place the device on a cold surface or immerse in a cold fluid.

**L. MATERIAL SPECIFICATIONS**  
Refer to the package label for the materials. This device is made of titanium or stainless steel.

**M. STORAGE CONDITIONS**  
Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

### N. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

If **non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.** Disinfection solutions should always be mixed to the manufacturer’s specification for concentration.

- After disinfection, the devices should be rinsed with cold distilled water or deionized sterile water. (Acceptable endotoxin levels for deionized water per AMI TR 34:2014 is < 10 EU/mL.)
- Dry devices thoroughly using compressed air, wipes, or an oven.
- Check devices for visible soil. Repeat cleaning if soil is visible and re-inspect.

### II. MACHINE (AUTOMATIC) CLEANING, AND THERMAL DISINFECTION

The instructions for Machine (Automatic) Cleaning and Thermal Disinfection may be followed as an alternative cleaning method to Manual Cleaning and Disinfection.

- Load the devices in the washer such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain.
- Run the automatic wash cycle. The minimum cycle parameters are listed below.

MINIMUM WASHING CYCLE PARAMETERS			
Phase	Recirculation Time	Temperature	Detergent
<b>Cold Pre-Wash</b>	3 Minutes	68 ± 9°F (20 ± 5°C)	NA
<b>Cleaning Wash</b>	10 Minutes	150 ± 9°F (65.5 ± 5°C)	Enzymatic or alkaline agent
<b>Rinse 1</b>	3 Minutes	122 ± 9°F (50 ± 5°C)	NA
<b>Rinse 2</b>	3 Minutes	122 ± 9°F (50 ± 5°C)	NA
<b>Thermal Disinfection Rinse</b>	5 Minutes	194 ± 9°F (90 ± 5°C)	NA
<b>Drying</b>	6 Minutes	239°F (115°C)	NA

- Automatic wash cleaning solutions can include, but are not limited to: needisher® Mediclean forte, Thermosep® alk clean, Prolysta® Ultra Concentrate Enzymatic Cleaner, and Proklean NpH Neutral Detergent. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™.**

If **non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**

- Check devices for visible soil. Repeat cleaning if soil is visible and re-inspect.

### K. STERILIZATION

This device may be provided either sterile or non-sterile. Check the package labeling for more information. For the devices that are not provided in a terminally sterilized configuration, sterilization is to be performed following cleaning, disinfection prior to use.

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