English

A. DEVICE DESCRIPTION

The Arthrex Low Profile Screws are headed and self-tapping. They are available as fully or partially threaded, and solid or cannulated. The screw filament ranges from 2.0 mm to 6.7 mm in diameter and from 8 mm to 120 mm in length (in 1, 2 or 5 mm increments).

B. INDICATIONS

The **Arthrex Low Profile Screws (2.0-3.0 mm solid)** are intended to be used as stand-alone bone screws, or in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, and hand. When used with a plate, the screw may be used with the Arthrex Low Profile, Small Fragment Plates, Distal Extremity Plates, Mesh Plates, and Distal Radius Plates.

The **Arthrex Low Profile Screws (2.0-3.0 mm cannulated)** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, and wrist.

The **Arthrex Low Profile Screws (3.5 mm and larger, solid)** are intended to be used as stand-alone bone screws, or in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur, and fibula. When used with a plate, the screws may be used with the Arthrex Low Profile, Small Fragment Plates, Fracture Plates, Distal Extremity Plates, Distal Radius Plates, Humeral Fracture Plates, and Osteotomy Plates.

The **Arthrex Low Profile Screws (3.5 mm and larger, cannulated)** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used with a plate, the screws may be used with the Arthrex Low Profile, Small Fragment Plates, Fracture Plates, Distal Extremity Plates, Distal Radius Plates, Humeral Fracture Plates, and Osteotomy Plates.

The **Arthrex Verted Cannulated Locking Screws (4.0 mm cannulated)** are intended to be used in a plate-screw system for internal bone fixation for bone fractures in the humerus. The screws are used with the Arthrex Humeral Fracture Plates. These screws are for EU distribution only.

The **Arthrex Blunt Tip Screws** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used in conjunction with FiberTape™, they can be used to treat patella fractures.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity: Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

E. WARNINGS

- An internal fixation device must never be reused.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be considered as temporary and should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- Devices that have been implanted for a long period of time may require the use of screw removal instrumentation. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of the intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult. (5) Pain, discomfort, or abnormal sensations due to the presence of the device. (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

F. MRI SAFETY INFORMATION

1. MR CONDITIONAL

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Low Profile Screws are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less
- Maximum MR system reported, with major averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Arthrex Low Profile Screws are expected to produce a maximum temperature rise of up to 14.9 °C after 15-minutes of continuous scanning
- Artifact Information**
In non-clinical testing, the image artifact caused by the Arthrex Low Profile Screws can extend up to approximately 120 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Use the appropriately sized drill bit for the screw.
- Damage to the driver or screw may result from failure to seat the driver fully into the screw or to align the driver properly with the screw.
- QuickFix™ Screws** Only: Bending the insertion post to remove it from the screw head is not recommended. Screws should be inserted by hand and not with powered equipment.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labels arrive intact.
- Contact Customer Service if the package has been opened or altered.

- All of the symbols used on the labeling along with the title, description and standard dimensions will either may be found on our website at www.arthrex.com/symbols/glossary.

I. VALIDATION

The recommended cleaning, disinfection, and sterilization methods in this DFU have been validated in compliance with federal and international guidance/standards. In accordance with ISO 17665, the half cycle “overall” approach was used for sterilization validation, and an approved a sterility assurance level (SAL) of 10⁻⁶. Cleaning, disinfecting, and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility/end user to use equipment cleared by the FDA or the local governing body and perform operations in accordance with directions from this DFU. In accordance with EN ISO 17664 during the validation process, limit values and a means for monitoring chemical residues following cleaning have been established for the product and have been analyzed utilizing cytotoxicity testing as the primary means of assessing safety of the detergent and/or disinfectant levels after the cleaning and washing processes have been carried out. The quality of water utilized for the cleaning and disinfection processes should be validated at the end-user site to ensure that residuals will not interfere with subsequent processing steps in accordance with AMI TR30-2011. Deionized water was the water quality used for manual and automated cleaning processing instructions performed to establish the required processing parameters as outlined in the “Cleaning” section of this document. The instructions in this DFU were developed using the guidance given in the following standards:

- ANSI/AMI S179, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.”
- ISO 17664: Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.
- ISO 17665-1: Sterilization of health care products – Part 1: Requirements for the sterilization process, including validation, and routine control of a sterilization process for medical devices
- ISO 10993-5: Biological evaluation of medical devices – Part 5: Tests for in vivo cytotoxicity

J. CLEANING

Certain Arthrex devices that may be used during this procedure are provided, non-sterile and must be adequately cleaned and sterilized prior to use re-use.

I. MANUAL CLEANING AND DISINFECTION

- The instructions for Manual Cleaning and Disinfection may be found in an alternative cleaning method to Machine (Automatic) Cleaning and Thermal Disinfection.
 - Immerse the device in an enzymatic or alkaline cleaning detergent solution. Cleaning solutions can include, but are not limited to: ENZO® enzymatic, neodisher® Mediclean forte, and Thermostat® alk clean. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylenepropylene), ABS (Acrylonitrilbutadiene Styrene), Utem™, Lexan™, and Cycloac™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**

L. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. This device is made of titanium or stainless steel.

M. STORAGE CONDITIONS

- Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date. Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.
- N. INFORMATION**
Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

II. MACHINE (AUTOMATIC) CLEANING, AND THERMAL DISINFECTION

The instructions for Machine (Automatic) Cleaning and Thermal Disinfection may be followed as an alternative cleaning method to Manual Cleaning and Disinfection.

- Load the devices in the washer such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain.
- Run the automatic wash cycle. The minimum cycle parameters are listed below.

MINIMUM WASHING CYCLE PARAMETERS				
Phase	Reirculation Time	Temperature	Detergent	
Cold Pre-Wash	3 Minutes	68 ± 9°F (20 ± 5°C)	N/A	
Cleaning Wash	10 Minutes	150 ± 9°F (65.5 ± 5°C)	Enzymatic or alkaline agent	
Rinse 1	3 Minutes	122 ± 9°F (50 ± 5°C)	N/A	
Rinse 2	3 Minutes	122 ± 9°F (50 ± 5°C)	N/A	
Thermal Disinfection Rinse	5 Minutes	194 ± 9°F (90 ± 5°C)	N/A	
Drying	6 Minutes	239°F (115°C)	N/A	

- Automatic wash cleaning solutions can include, but are not limited to: neodisher® Mediclean forte, Thermostat® alk clean, Physiolys® Ultra Concentrate Enzymatic Cleaner, and ProKlean NpI Neutral Detergent. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylenepropylene), ABS (Acrylonitrilbutadiene Styrene), Utem™, Lexan™, and Cycloac™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**
- Check devices for visible soil. Repeat cleaning if soil is visible and re-inspect.

K. DECONTAMINATION

This device is provided sterile or non-sterile. Check the package labeling for more information. For the devices that are not provided a terminally sterilized configuration, sterilization is to be performed following cleaning, disinfection prior to use. This device (i.e. for open but unused) may be re-sterilized following cleaning, disinfection prior to use.

Sterilization wrap that has been FEA-cleared for the specified sterilization cycles should be used. Ensure that the wrap is large enough to contain the device without stressing the seals.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to re-use. Please refer to the Arthrex Instruments Directions for Use and ANSI/AMI S179, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities” for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:				
	Exposure Temperature	Exposure Time	Drying Time	
	121°C (250°F)	30 Minutes	15 to 30 Minutes	
Gravity-Displacement Steam Sterilization Cycle	132°C (270°F)	15 Minutes	15 to 30 Minutes	
	135°C (275°F)	10 Minutes	30 Minutes	
Pre-vacuum Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes	
	135°C (275°F)	3 Minutes	16 Minutes	

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:

	Exposure Temperature	Exposure Time	Drying Time	
	132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes	
Gravity-Displacement Steam Sterilization Cycle	121°C (250°F)	30 Minutes	15 to 30 Minutes	
Pre-vacuum Cycle	132°C – 135°C (270°F – 275°F)	4 Minutes	30 to 30 Minutes	

Cooling – The device must be adequately cooled, after being removed from the sterilizer. It should not be touched during the cooling process. Do not place the device on a cold surface or immerse in a cold fluid.

L. MATERIAL SPECIFICATIONS
Refer to the package label for the materials. This device is made of titanium or stainless steel.

M. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date. Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

N. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

3. All of the symbols used on the labeling along with the title, description and standard dimensions will either may be found on our website at www.arthrex.com/symbols/glossary.

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