### Symbols Used on Labeling

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## English DFU-0274-1r0

### A. DEVICE DESCRIPTION

The Angel® PRP (Platelet Rich Plasma) Kit contains blood draw components inclusive of an infusion set, tourniquet, alcohol pad, gauze, Band-Aid, 60mL syringes, in addition to the Angel cPRP (concentrated Platelet Rich Plasma) Processing Set.

### B. INDICATIONS

The Angel PRP Kit is intended to facilitate blood draws for use with the Angel cPRP Processing Set in the production of PRP.

The Angel cPRP Processing Set is intended for use with the Angel Concentrated Platelet Rich Plasma (cPRP) System to separate and collect an autologous plasma fraction rich in platelets and red blood cells from the patient's whole blood or a small mixture of blood and bone marrow perioperative to a surgical procedure.

### C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.
3. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
4. Do not use for surgeries other than those indicated.

### D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Foreign body reactions.
3. Hematoma.

### E. WARNINGS

1. Do not re-sterilize this device.
2. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
3. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.
4. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.

### F. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also

provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

### G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.

### H. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

### I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The contents are manufactured from ABS (acrylonitrile butadiene styrene), Polycarbonate, Polyethylene, Polypropylene, Polyurethane, PVC (polyvinyl chloride), Stainless Steel, and Isoprene Rubber.

### J. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

### K. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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