# English

**In the USA**

Angel® cPRP System with Aspiration Kit

**A. DEVICE DESCRIPTION**

The Angel® cPRP System with Aspiration Kit contains a fenestrated bone marrow aspiration needle, syringe, clot filter, needle, cup, and an Angel Concentrated Platelet Rich Plasma (cPRP) Processing Set.

The devices comprised of the Angel cPRP System with Aspiration Kits are designed to be used with the Angel Concentrated Platelet Rich Plasma (cPRP) System to deliver PRP concentrate from bone marrow aspirate (BMA).

**B. INDICATIONS**

**For the USA:** The Angel Concentrated Platelet Rich Plasma (cPRP) System is intended to be used in the clinical laboratory or intraoperatively at the point-of-care for the safe and rapid preparation of platelet poor plasma and concentrate (platelet rich plasma) from a small sample of whole blood or a small mixture of blood and bone marrow. The platelet rich plasma from the Angel Concentrated Platelet Rich Plasma (cPRP) System can also be mixed with autograf and/or allograft bone prior to application to an orthopedic site.

**For outside the USA:** The devices comprised of the Angel cPRP & Bone Marrow Aspirate (BMA) Kits are intended to be used with the Angel Concentrated Platelet Rich Plasma (cPRP) System to deliver customized bone marrow concentrate from a small mixture of blood and bone marrow for the treatment of osseous defects within the hip/ acetabulum.

**C. CONTRAINDICATIONS**

- Blood supply limitations and previous infections, which may retard healing.
- Any active infection or blood supply limitations.
- Conditions which tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- Do not use for surgeries other than those indicated.

**D. ADVERSE EFFECTS**

- Infections, both deep and superficial.
- Foreign body reactions.
- Hematoma.

**E. WARNINGS**

- Do not re-sterilize this device or any of its components.
- Do not reuse this device or any of its components.
- Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
- A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

**F. PRECAUTIONS**

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

**G. PACKAGING AND LABELING**

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

**H. STERILIZATION**

The device is provided sterile. Refer to the package label for the sterilization method.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities”

**I. MATERIAL SPECIFICATIONS**

Refer to the package label for the materials.

- Needle:** Surgical-grade stainless steel.
- Syringe:** Polypropylene with a polyisoprene seal.
- Preparation Tray:** Polypropylene.

**J. STORAGE CONDITIONS**

This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

**K. INFORMATION**

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

**L. DIRECTIONS FOR USE**

Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require more comprehensive aspiration procedures.

**I. BASIC BONE MARROW ASPIRATION PROCEDURE**

- Locate and prepare the bone marrow aspiration entry site.

- Insert the aspiration needle.
- Attach the male connector of the syringe to the female connector handle of the aspiration needle.
- Aspirate.
- Transfer from the syringe as deemed necessary.

**Outside USA**

Angel® cPRP &amp; Bone Marrow Aspirate (BMA) Kit

**A. DEVICE DESCRIPTION**

The Angel® cPRP & Bone Marrow Aspirate (BMA) Kit contains a fenestrated bone marrow aspirate needle, syringe, clot filter, needle, cup, and Angel Concentrated Platelet Rich Plasma (cPRP) Processing Set..

The devices comprised of the Angel cPRP & Bone Marrow Aspirate (BMA) Kits are designed to be used with the Angel Concentrated Platelet Rich Plasma (cPRP) System to deliver PRP concentrate from bone marrow aspirate (BMA).

**B. INDICATIONS**

**For the USA:** The Angel Concentrated Platelet Rich Plasma (cPRP) System is intended to be used in the clinical laboratory or intraoperatively at the point-of-care for the safe and rapid preparation of platelet poor plasma and concentrate (platelet rich plasma) from a small sample of whole blood or a small mixture of blood and bone marrow. The platelet rich plasma from the Angel Concentrated Platelet Rich Plasma (cPRP) System can also be mixed with autograf and/or allograft bone prior to application to an orthopedic site.

**For outside the USA:** The devices comprised of the Angel cPRP & Bone Marrow Aspirate (BMA) Kits are intended to be used with the Angel Concentrated Platelet Rich Plasma (cPRP) System to deliver customized bone marrow concentrate from a small mixture of blood and bone marrow for the treatment of osseous defects within the hip/acetabulum.

**C. CONTRAINDICATIONS**

- Blood supply limitations and previous infections, which may retard healing.
- Any active infection or blood supply limitations.
- Conditions which tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- Do not use for surgeries other than those indicated.

**D. ADVERSE EFFECTS**

- Infections, both deep and superficial.
- Foreign body reactions.
- Hematoma.

**E. WARNINGS**

- Do not re-sterilize this device or any of its components.
- Do not reuse this device or any of its components.
- Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
- A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

**F. PRECAUTIONS**

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration

**G. PACKAGING AND LABELING**

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

**H. STERILIZATION**

The device is provided sterile. Refer to the package label for the sterilization method.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities”

**I. MATERIAL SPECIFICATIONS**

Refer to the package label for the materials.

- Needle:** Surgical-grade stainless steel.
- Syringe:** Polypropylene with a polyisoprene seal.
- Preparation Tray:** Polypropylene.

**J. STORAGE CONDITIONS**

This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

**K. INFORMATION**

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

**L. DIRECTIONS FOR USE**

Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require more comprehensive aspiration procedures.

**I. BASIC BONE MARROW ASPIRATION PROCEDURE**

- Locate and prepare the bone marrow aspiration entry site.
- Insert the aspiration needle.
- Attach the male connector of the syringe to the female connector handle of the aspiration needle.
- Aspirate.
- Transfer from the syringe as deemed necessary.

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