English

A. DEVICE DESCRIPTION

Plates are available in different shapes, sizes and orientations (e.g. left and right types). The plates have specified sized holes for screws to provide fixation. Some plates with multiple holes (i.e. Arthrex Mesh Plates) may be cut to size.

The Distal Clavicle Plate Button is designed to fit securely in the holes of the fracture plates. The Plate Button is available in a size only. The Distal Clavicle Plate Button is available as a single button or assembled in a construct with Arthrex suture washers. Washers are available in diameters of 6 – 13 mm.

B. INDICATIONS

The *Arthrex Fracture Plates and Ankle Fusion Plates* are intended to be used for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus and fibula. The *Arthrex Fracture Plates* are to be used with the *Arthrex Low Profile Screws*.

The *Distal Clavicle Plate Button* is intended for use with the clavicle plates for clavicle indications such as for the treatment of syndesmotic trauma, such as fixation of acromioclavicular separations due to coracoclavicular ligament disruption, and this button may not be used alone. The button is intended to be used with #5 FiberWire® or FiberTap®.

The *Arthrex Low Profile Plates, Arthrex Compression Plates and Calcaneal Fracture Plates* are intended to be used for internal bone fixation for bone fractures, fusions, or osteotomies in the ankle, foot, hand, and wrist, such as opening wedge osteotomies of Hallux Valgus.

The *Arthrex Distal Extremity Plate System and Arthrex Mesh Plates* are intended for use in stabilization of fresh fractures, revision procedures, osteotomies, joint fusion and reconstruction of small bones and bone fragments of the hand/wrist, ankle/ankle, and osteopic bone.

The *Arthrex Distal Radius System* is intended for fixation of fractures, fusions, osteotomies and non-unions of the radius, ulna, olecranon, metacarpal, metatarsal and malloles.

The *Arthrex* suture washers and certain plates are intended as adjunct devices.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infections or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this metal device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

E. INTERNAL FIXATION DEVICES MAY NEVER BE REUSED

- All metallic implant devices must never be reused.
- All metallic implant devices used for this surgical procedure should always have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection of this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

F. MRI SAFETY INFORMATION

- This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Do not bend the plate near the locking hole. Bending the plate near the locking hole can distort the holes through, which prohibits insertion of the screw.
- Repeated bending of the plate at the same location, or by creating excessive acute angles may potentially lead to premature plate fatigue, failure and or breakage in situ.
- Screws should be inserted by hand and not with powered equipment.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbology.aspx.

I. VALIDATION

The recommended cleaning, disinfection, and sterilization methods in this DPU have been validated in compliance with federal and international guidance standards. In accordance with ISO 17665, the half cycle "overkill" approach was used for sterilization validation, and demonstrates a sterility assurance level (SAL) of 10⁻⁶. Cleaning, disinfecting, and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility end-user to use equipment cleared by the FDA or the local governing body and perform operations in accordance with directions from this DPU.

In accordance with EN ISO 17664 during the validation process, limit values and a minimum decontamination level (SAL) following cleaning have been established for the product and have been analyzed utilizing cytotoxicity testing as the primary means

of assessing safety of the detergent and/or disinfectant levels after the cleaning and washing processes have been carried out.

The quality of water utilized for the cleaning and disinfection processes should be validated at the end-user site to ensure that residuals will not interfere with subsequent processing steps in accordance with AAMI TR30:2011. Deionized water was the water quality used for manual and automated cleaning process

validations performed to establish the required processing parameters as outlined in the "Cleaning" section of this document. The instructions in this DPU were developed using the guidance given in the following standards:

- ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities"
- ISO 17664: Sterilization of medical devices – Information to be provided by the manufacturer for the processing of reutilizable medical devices
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices
- ISO 10993-5: Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity

J. CLEANING

Certain Arthrex devices that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use.

I. MANUAL CLEANING AND DISINFECTING

The instructions for Manual Cleaning and Disinfection may be followed as an alternative cleaning method to Machine (Automatic) Cleaning and Thermal disinfection.

- Immerse the device in an enzymatic or alkaline cleaning detergent solution. Cleaning solutions can include, but are not limited to: ENZOL® enzymatic, needisher® Mediclean forte, and Thermosep® alkali clean. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.** Cleaning solutions should always be mixed to the manufacturer's specification for concentration and cleaning shall be conducted at ambient temperature unless otherwise stated in the cleaning solution manufacturer's instructions.

- Scrub device with a soft brush, paying special attention to areas where debris might collect. Solutions can include, but are not limited to, CIDEX®, WAVIDE®-01, Gigasep®, Kohrsoln®, and equivalent products. Use the supplier's instructions for preparing the solution. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™. If non-neutral pH disinfection chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.** Disinfection solutions should always be mixed to the manufacturer's specification for concentration.
- After disinfection, the devices should be rinsed with cold distilled water or deionized sterile water. (Acceptable endotoxin levels for deionized water per AAMI TIR 34:2014 is <10 EU/ml.)
- Dry devices thoroughly utilizing compressed air, wips, or cloths.
- Check devices for visible soil. Repeat cleaning if soil is visible and re-insept.

- Immerse devices in disinfection solutions for a minimum of 20 minutes. Suitable disinfection solutions can include, but are not limited to: CIDEX®, WAVIDE®-01, Gigasep®, Kohrsoln®, and equivalent products. Use the supplier's instructions for preparing the solution. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™. If non-neutral pH disinfection chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.** Disinfection solutions should always be mixed to the manufacturer's specification for concentration.

- After disinfection, the devices should be rinsed with cold distilled water or deionized sterile water. (Acceptable endotoxin levels for deionized water per AAMI TIR 34:2014 is <10 EU/ml.)
- Dry devices thoroughly utilizing compressed air, wips, or cloths.
- Check devices for visible soil. Repeat cleaning if soil is visible and re-insept.

- Automatic wash cleaning solutions can include, but are not limited to: needisher® Mediclean forte, Thermosep® alkali clean, Prolystica® Ultra Concentrate Enzymatic Cleaner, and ProkZen NpH Neutral Detergent. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**

- Check devices for visible soil. Repeat cleaning if soil is visible and re-insept.
- Run the automatic wash cycle. The minimum cycle parameters are listed below:

MINIMUM WASHING CYCLE PARAMETERS			
Phase	Recirculation Time	Temperature	Detergent
Cold Pre-Wash	3 Minutes	68 ± 9°F (20 ± 5°C)	N/A
Cleaning Wash	10 Minutes	150 ± 9°F (65.5 ± 5°C)	Enzymatic or alkaline agent
Rinse 1	3 Minutes	122 ± 9°F (50 ± 5°C)	N/A
Rinse 2	3 Minutes	122 ± 9°F (50 ± 5°C)	N/A
Thermal Disinfection Rise	5 Minutes	194 ± 9°F (90 ± 5°C)	N/A
Drying	6 Minutes	239°F (115°C)	N/A

- Automatic wash cleaning solutions can include, but are not limited to: needisher® Mediclean forte, Thermosep® alkali clean, Prolystica® Ultra Concentrate Enzymatic Cleaner, and ProkZen NpH Neutral Detergent. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**

K. STERILIZATION
This device may be provided either sterile or non-sterile. Check the package labeling for more information. For the devices that are not provided in a terminalized sterilization configuration, sterilization is to be performed following cleaning, disinfection prior to use. This device (i.e. for open but unused) may be re-sterilized following cleaning, disinfection prior to use.

- Check devices for visible soil. Repeat cleaning if soil is visible and re-insept.
- Automatic wash cleaning solutions can include, but are not limited to: needisher® Mediclean forte, Thermosep® alkali clean, Prolystica® Ultra Concentrate Enzymatic Cleaner, and ProkZen NpH Neutral Detergent. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**

K. STERILIZATION
This device may be provided either sterile or non-sterile. Check the package labeling for more information. For the devices that are not provided in a terminalized sterilization configuration, sterilization is to be performed following cleaning, disinfection prior to use. This device (i.e. for open but unused) may be re-sterilized following cleaning, disinfection prior to use. Sterilization wrap that has been FDA cleared for the specified sterilization cycles should be used. Ensure that the wrap is labeled enough to contain the device without stressing the seals.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to the Arthrex Instruments Directions for Use and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

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