English

A. DEVICE DESCRIPTION

The PASTA Bridge Repair Implant System consists of multiple biocomposite or PEEK anchors and manual surgical instruments for use in PASTA Bridge Repair.

B. INDICATIONS

The SutureTak®, Corkscrew and SwiveLock® Suture Anchors are intended to be used for suture or tissue fixation in the shoulder.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infection, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be reused.
- Do not re-sterilize this device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of the device, the patient leaflet and the patient implant card should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
- Attempting implantation into hard, dense bone and/or drilling/punching smaller diameter holes than

recommended may cause failure (breakage) of the implant during insertion.

- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Safe

The Corkscrew, SwiveLock, and SutureTak suture anchors devices manufactured from only polyester, polyetheretherketone (PEEK) poly (L-Lactide, PLLA), Poly (L-Lactide-co-D acid, PLDLA), and/or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP) are MR safe.

G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate suture-anchor type and size based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Insert the anchor with the same orientation as that of the prepared drilled bone hole to avoid damage to the anchor.
- Make sure to use the recommended drill bit or punch to create the bone socket.
- SwiveLock suture anchors only:** Ensure that the anchor body is in full contact with the bone before advancing the anchor body into the prepared bone socket.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Anchor: Biocomposite anchors are made of Poly (L-Lactide) or Poly (L-Lactide-co-D, L-Lactide) and tricalcium phosphate (TCP).

Some anchors consist of either one or two components. Each component is manufactured from either/ or polyetheretherketone (PEEK), Poly (L-Lactide acid, PLLA), Poly (L-Lactide-co-D acid, PLDLA), or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP).

Suture: See package label for size and type of suture provided with device.

FiberWire®, FiberWire CL, TigerWire®, TigerWire CL, FiberTape®, and TigerTape™ sutures are made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. Additional materials may include silicone elastomer coating (except suture with the suffix -Tape), cyanoacrylate, and may include nylon.

SutureTape suture is made from Ultra High Molecular Weight Polyethylene (UHMWPE), and polyester. Additional materials may include nylon and/or cyanoacrylate. The #2 Suture and #5 Suture are a polyblend made from Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester.

The sutures supplied with the suture anchor meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The suture dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon.

K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C/ 90° F, and should not be used after the expiration date.

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