Symbols Used on Labeling
Auf der Verpackung verwendete Symbole
Symboles utilisés sur l'étiquette
Simboli usati sull'etichetta
Símbolos utilizados en las etiquetas
Símbolos utilizados na rotulagem

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English

DFU-0054-2r0

A. DEVICE DESCRIPTION
These devices consist of the Bio-FASTak®, FASTak™, SutureTak™ and FiberTak™ Suture Anchors.

These suture anchor families consist of rigid and threaded implants and all-suture implants (i.e. FiberTak). FiberTak suture anchors are comprised of a polyester sheath with single or multiple sutures assembled through the sheath. All implants are provided pre-loaded on an inserter and may include eyelets which are preloaded with various types of suture, with or without needles.

B. INDICATIONS
The Bio-FASTak, FASTak, SutureTak and FiberTak Suture Anchors are intended to be used for suture or tissue fixation in the foot, ankle, knee, hand, wrist, elbow, shoulder, and hip (devices with FiberWire® or #2 polyblend Suture).

Surgeons must apply their professional judgment when determining the appropriate suture-anchor type and size based on the specific indication, preferred surgical technique, and patient history.

The Bio-FASTak, FASTak, SutureTak and FiberTak Suture Anchors are intended for the following indications:

Shoulder: Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulolabral Reconstruction
Foot/Ankle: Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Mid-foot Reconstruction, Metatarsal Ligament Repair, (**SutureTak and FiberTak only**): Digital Tendon Transfers.

Knee: Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis

Hand/Wrist: Scapholunate Ligament Reconstruction, Repair/Reconstruction of Collateral Ligaments, (**SutureTak and FiberTak only**): Carpal Ligament Reconstruction, Repair of Flexor and Extensor Tendons at the PIP, DIP and MCP joints for all digits, Digital Tendon Transfers. (**Knotted FiberTak and Micro BioComposite SutureTak only**): Carpometacarpal joint arthroplasty (basal thumb joint arthroplasty).

Elbow: Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction
Hip (devices with FiberWire® or #2 Suture {polyblend only}): acetabular labral repair, Capsular repair (**excludes Knotless FiberTak**). Exclusions to hip indications include knotted FiberTak suture anchors, and 2.4 mm FASTak.

The 2.0mm and 2.4mm Hip SutureTaks are indicated for acetabular labral repair **ONLY**.

C. CONTRAINDICATIONS

1. Insufficient quantity or quality of bone.
2. Blood supply limitations and previous infections, which may retard healing.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
4. Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
5. Any active infection or blood supply limitations.
6. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
7. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
8. Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

1. Infection, both deep and superficial.
2. Foreign body reactions.
3. Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
4. **FASTak only:** Shoulder dislocation/subluxation.

E. WARNINGS

1. An internal fixation device must never be reused.
2. Do not resterilize this device.
3. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
4. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
5. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
6. Detailed instructions on the use and limitations of the device should be given to the patient.
7. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
8. Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

F. MRI SAFETY INFORMATION

1. This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.

G. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides

detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

2. **FASTak only:** Stopping and restarting the drill may result in excess torque being applied to the implant, which may cause the device to fail.

3. **FASTak only:** It is important to stop drilling as soon as the chuck contacts the guide or grasper. Failure to do so may result in damage to the suture and/or implant.

4. Insert the anchor with the same orientation as that of the prepared drilled bone hole to avoid damage to the anchor.

5. **SutureTak and FiberTak only:** Anchors Loaded on Flexible Drivers: The drill guide tip must remain in contact with the bone surface during the drilling and implant impaction steps of the procedure. Failure to do so may result in difficulty seating the implant to its intended depth.

H. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.

I. STERILIZATION
This device is provided sterile. Refer to the package labeling for the sterilization method.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

J. MATERIAL SPECIFICATIONS
Refer to the package label for the materials.

Anchor: The device is made of titanium, PEEK (polyetheretherketone), Poly (L-Lactide) or Poly (L-Lactide-co-D, L-Lactide). FiberTak is made with a polyester sheath.

Biocomposite anchors are made of Poly (L-Lactide) or Poly (L-Lactide-co-D, L-Lactide) and tricalcium phosphate (TCP).

Suture: See package label for size and type of suture provided with device.

FiberWire, FiberWire CL, TigerWire, TigerWire CL, FiberTape, and TigerTape sutures are made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. Additional materials may include silicone elastomer coating (except suture with the suffix -Tape), cyanoacrylate, and may include nylon.

LabralTapeTM is made from Ultra High Molecular Weight Polyethylene (UHMWPE) and may include polyester or nylon.

SutureTape™ is made from Ultra High Molecular Weight Polyethylene (UHMWPE), and polyester. Additional materials may include nylon and/or cyanoacrylate.

The #2 suture is a polyblend made from Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester.

or
The suture is made from Polyester, polytetrafluoroethylene (PTFE) coating, cyanoacrylate.

The sutures supplied with the suture anchor meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The suture dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon.

K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-Bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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