<p><b>Symbols Used on Labeling</b></p> <p>Auf der Verpackung verwendete Symbole</p> <p>Symboles utilisés sur l'étiquetage</p> <p>Simboli usati sull'etichetta</p> <p>Símbolos utilizados en las etiquetas</p> <p>Símbolos utilizados na embalagem</p>
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

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## English

### A. INDICATIONS

The Chondral Dart is intended for the use in the fixation of small bone fragments, such as apical fragments, osteochondral fragments, and cancellous fragments in the knee. Specific applications include the following: Apical fragments (radial head, patellar rim, navicular, metacarpal/metatarsal), osteochondral fragments (talus vault, femoral condyle) and cancellous fragments (talus).

The Meniscal Dart is intended for the repair of meniscal tears that would otherwise be considered for standard repair using suture.

### B. CONTRAINDICATIONS

- Chondral Dart only:** Insufficient quantity or quality of bone.
- Meniscal Dart only:** Surgeries other than those indicated such as flap tears, or complex bucket handle tears.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb, or disrupt the growth plate.

### C. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

### D. WARNINGS

- Meniscal Dart only:** Meniscal Darts may protrude through the joint capsule. It is important to use a Dart of the appropriate length. Protruding portions of the Dart may be cut off through small separate incisions, if necessary.
- An internal fixation device must never be reused.
- Do not re-sterilize this device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating or migration in the MR environment. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

### E. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Meniscal Dart only:** Do not apply excessive posterior pressure to the sheath.
- Meniscal Dart only:** If necessary, cut off any dart length protruding through the joint capsule.

### F. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.

### G. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

## DFU-0052r I I

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities", for specific information.

### H. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.  
The Chondral Dart is made of Poly L-lactide.  
The Meniscal Dart is made of Poly (L-lactide-co-D, L-lactide).

### I. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at a maximum 32° C/90° F, and should not be used after the expiration date.

### J. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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