English

A. DEVICE DESCRIPTION

The FibuLock® Nail is designed to provide fixation of fibula fractures while they heal.

B. INDICATIONS

The FibuLock System is intended for use in the fixation of fibula fractures and osteotomies.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone/soft tissue.
- Patients with poor prognosis for healing.
- Blood supply limitations and previous infections, which may retard healing.
- Nickel sensitivity. If suspected, tests are performed prior to implantation.
- Active or latent infection.
- Patients who are unwilling or incapable of following postoperative care instructions.
- Skeletal immaturity if the surgical technique would affect the development of the bone after placement of the fixation implant.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Potential adverse events related to any surgery: Reactions to anesthesia, the anesthetic or other medications, bleeding, infection (both deep and superficial), ileus, blood vessel damage, nerve or soft tissue damage, atelectasis, pneumonia, hematoma, seroma, wound dehiscence or incisional hernia, urologic problems, embolism, anemia, colitis, thrombophlebitis, heart attack, stroke or death.
- Potential adverse events related to the device: Bending, cracking, or fracture of the implant, failure to achieve fracture healing, or loosening of the implant. Limb shortening or loss of anatomic position with nonunion, malunion with rotation or angulation. Irritationl injury of soft tissues, including impingement syndrome, or tissue reactions which include macrophage and foreign body reactions adjacent to implants, Metal sensitivity reactions and/or allergic reactions to foreign materials, pain, discomfort, or abnormal sensations due to the presence of the implant. Necrosis of bone or bone resorption. Necrosis of the tissue or inadequate healing. Generation of material debris possibly resulting in bodily response. Additional surgery may be necessary for implant removal, repositioning or replacement.

E. WARNINGS

- Do not re-sterilize this device.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Bone fixation devices are neither intended to carry the full load of the patient, nor intended to carry a significant portion of the load for extended periods of time. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration.
- The correct selection and placement of the implant is extremely important. The appropriate type and size should be selected for the patient. Failure to use the correct implant size or improper positioning may result in loosening, bending, cracking, or fracture of the device, or bone or both.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- The patient must be cautioned about the use, limitations, and possible adverse effects of this implant including the possibility of the device failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight and/or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing.

F. MRI SAFETY INFORMATION

- This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.

G. PRECAUTIONS

- Implants are single-use only and should never be reused.
- Use care in handling and storage of implant components. Cutting, sharply bending or scratching the surface can significantly reduce the strength and fatigue resistance of the implant system.

- The surgeon should be thoroughly familiar with the devices, instruments and surgical technique before surgery. A detailed surgical technique is available for surgeon reference. Proper surgical procedures and techniques are the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of the procedures based on his or her personal medical training and experience. During the post-operative phase, in addition to mobility and muscle training, it is of particular importance that

the physician keeps the patient well informed about post-surgical activity requirements.

- Implants should not be forced, twisted, bent, or pried during insertion as this may damage, break, or bend the implant and render it unusable. Damaged implants should not be implanted.
- The use of locking screws is necessary for fracture stability. Refer to the surgical technique or product listing for the available screw sizes.
- The surgeon should take care with patient selection when treating patients with preexisting conditions that may affect the success of the surgical procedure.
- The patient should be advised that an additional procedure for the removal of the implant is sometimes necessary. Providing postoperative instructions to patients and appropriate nursing care is critical. Early weight and/or load bearing substantially increases implant loading and increases the risk of loosening, bending or breaking the device. Patients who are obese and/or noncompliant, as well as patients who could be pre-disposed to delayed union or non-union, must have auxiliary support. Patients should be cautioned against unassisted weight-bearing activity without physician direction or medical release. Postoperative care and physical therapy should be structured to prevent loading of the operative extremity until directed by the physician.
- Even after full healing, the patient should be cautioned that re-fracture is possible and that a revision surgery may be necessary.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbolsglossary**.

I. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities” for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. The FibuLock Nail is manufactured from medical grade stainless steel.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

L. REVISION SURGERY / IMPLANT REMOVAL

Metallic implants can loosen, fracture, corrode, migrate, cause pain, or stress shield bone even after a fracture is healed. If any of these complications occur, the surgeon must make the final decision on implant removal.

M. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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