<p>Symbols Used on Labeling</p> <p>Auf der Verpackung verwendete Symbole</p> <p>Symboles utilisés sur l'étiquetage</p> <p>Simboli usati sull'etichetta</p> <p>Símbolos utilizados en las etiquetas</p> <p>Símbolos utilizados na embalagem</p>

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English	DFU-0162VGrI
A. DEVICE DESCRIPTION The ABPS device is comprised of a unique centrifuge tube that contains specially designed glass beads that maximize the surface-area-to-volume ratio inside the tube. After incubation, the autologous serum produced contains elevated levels of anti-inflammatory cytokines. The Arthrex ABPS kit contains the Arthrex ABPS device, a 60 mL syringe, and a winged blood recovery needle.	7. Place the Arthrex ABPS device horizontally on the middle rack in incubator for approximately 24 hours at 37°C, not to exceed 24 hours. 8. After incubation, let the device stand for 20 minutes then place the Arthrex ABPS device into a centrifuge for 10 minutes at 4000 rpm. Make sure to balance centrifuge with an appropriate counter weight. 9. Remove white tethered cap. Slowly draw the serum into a 20mL syringe using a spinal needle, being careful not to pull up red blood cells. 10. Attach a sterile 0.22µm filter between the 20mL syringe containing the serum and an empty, sterile 6mL syringe. Transfer 4mL of serum to the 6mL syringe, detach the 6mL syringe and cap. Two serum filtrations are recommended.
B. INDICATIONS The Arthrex ABPS is indicated for animals exhibiting pain from mild to moderate osteoarthritis. The serum is derived from the animal's own blood and reduces the possibility of adverse allergic or anaphylactic side effects.	11. Repeat until all the serum has been transferred through the sterile filter into the 6mL syringes. Individual doses may be used immediately or frozen at or less than -18°C for administration at a later date. 12. Excess blood and/or serum should be disposed of appropriately.
C. CONTRAINDICATIONS 1. Reuse. 2. Animals with severe allergies. 3. Use on infected sites. 4. Joint injection other than those described in the Indications.	

D. ADVERSE EFFECTS
1. Infection. 2. Fever, pain, local sensation, red flare, and inflammation. 3. Hematoma. 4. Cellulitis.

E. WARNINGS
1. The use of Arthrex ABPS must be decided based on thorough consideration of the animal's condition and the condition of the osteoarthritic area. 2. Arthrex ABPS is supplied sterile. After unsealing the package, the devices should be handled under sterile conditions and used immediately. 3. Do not overfill the Arthrex ABPS construct. 4. Arthrex ABPS is intended for use by licensed veterinarians 5. Do not microwave serum

F. PACKAGING AND LABELING
1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact. 2. Contact Customer Service if the package has been opened or altered.

G. STERILIZATION
The device is supplied sterile. Refer to the package label for the sterilization method.

H. MATERIAL SPECIFICATIONS
Refer to the package label for the materials. Arthrex ABPS Construct: Polypropylene, with tethered polyethylene male Luer caps 60 mL Syringe: Polypropylene Winged Blood Recovery Needle Set: <i>Needle</i> , stainless steel; <i>Protector and connector cap</i> , polyethylene; <i>Connector</i> , polymethylmethacrylate; <i>Wing and tubing</i> , TOTM PVC (Tris (2-Ethylhexyl) Trimellitate) polyvinyl chloride. Beads: Specially designed Glass

I. STORAGE CONDITIONS
This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

J. INFORMATION
For more information, or a demonstration, contact your local Arthrex representative.

K. DIRECTIONS FOR USE
Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.

- Using aseptic technique, prepare the jugular vein for needle penetration. Attach the enclosed butterfly needle to the 60mL syringe and puncture the jugular vein.
- Harvest ~50mL of whole blood into the 60mL syringe.
- Remove needle from jugular and detach from syringe.
- Unscrew the red tethered cap and loosen the white tethered cap on the Arthrex ABPS construct.
- Attach the syringe with blood to the red luer and inject blood into the Arthrex ABPS device, holding the syringe and Arthrex ABPS device at an angle of approximate ~60°. The blood should flow against the side of the construct. (Take care not to inject above *Max Fill Line*.)
- Remove syringe, recap and tighten both tethered male caps, and gently rock the blood in the Arthrex ABPS device, inverting at least 10 times.

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Arthrex Vet Systems, GmbH