<b>English</b>
<b>A. DEVICE DESCRIPTION</b>
The articular cobalt chrome or UHMWPE (Ultra High Molecular Weight Polyethylene) surface functions as the main joint surface replacement and provides a similar surface geometry as the removed cartilage. The titanium mesh on the back supports bone ingrowth. The threaded post is fixed in a drill hole within the cancellous bone and provides primary fixation for this device to promote bone ingrowth.

### B. INDICATIONS

The Arthrex Partial Eclipse is intended to be used for the repair of a painful and /or functionally insufficient shoulder joint resulting from degenerative disease due to primary (osteoarthritis), posttraumatic, post inflammatory, congenital or septic degeneration. Specifically, the Arthrex Partial Eclipse is indicated for use in the replacement of articular or full thickness cartilage (focal chondral) defects.

Caution: Federal law restricts this device to sale by or on the order of a physician.

### C. CONTRAINDICATIONS

- Non-localized defects.
- Non-continuous defects.
- Defects larger than available implant.
- Inflammatory degenerative joint disease.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Blood supply limitations and previous infections, which may retard healing.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb, or disrupt the growth plate.
- Pathological conditions, such as insufficient quantity or quality of bone (e.g., cystic changes or severe osteopenia), which may compromise implant fixation.
- Do not use for surgeries other than those indicated.

### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Loosening, migration or loss of fixation of implant.
- Fretting and crevice corrosion can occur at the interface between the implant components.
- Fatigue fracture of the implants as a result of bone resorption around the implant components.
- Wear and damage to the implant articulating surface.
- Intraoperative or postoperative bone fracture.
- Allergic-like reactions to materials.

### E. WARNINGS

- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- Proper patient selection is critical to outcomes. Use of the Partial Eclipse on patients with contraindications may impair implant function or lead to poor outcomes.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device. Patient sensitivity to the device materials should be considered prior to implantation.
- Detailed instructions on the use and limitations of this device, the patient leaflet and the patient implant card should be given to the patient.
- Do not resterilize this device.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating or migration in the MR environment. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.

- Follow your institutions policy for safe disposal of all needles and other sharps or medical waste.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

### F. PRECAUTIONS

Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
The treated bone should be of sufficient bone stock to support loading and show a conforming surface to the implant shape.

#### I. GUIDE PLACEMENT

- Take care at all times to identify and protect nerve structures
- Before placing the drill sleeve against the bone, soft tissue should be bluntly dissected down to bone with potential nerve structures being identified and protected at all times.
- Place the Soft Tissue Protector over the Guide Pin to avoid potential injury to potential nerve structures and other soft tissue.
- Take care when inserting the Reamer and during reaming to avoid contact with the corresponding joint surface.
- Mate the Retro-Pin to the Reamer by hand prior to connecting to a power system.
- Ream the articular surface until the Reamer’s depth stop fully engages.

#### III. IMPLANTATION

- Take care when inserting the Partial Eclipse or Partial Glenoid UHMWPE device to avoid contact with the corresponding joint surface.
- The Partial Eclipse or Partial Glenoid UHMWPE device must be fully seated into the previously reamed articular surface.

### G. MRI SAFETY INFORMATION

#### 1. MR Conditional

*Non-clinical testing demonstrated that the Partial Eclipse is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Partial Eclipse is expected to produce a maximum temperature rise of up to 4.6° C after 15-minutes of continuous scanning.

##### a. Artifact Information

*In non-clinical testing, the image artifact caused by the Partial Eclipse can extend up to approximately 16 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

### H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

### I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instrumnts that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

### J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

**Articulating Surface:** Cobalt Chrome or Ultra High Molecular Weight Polyethylene (UHMWPE)
**Fixation Surface:** CP (commercially pure) Titanium porous coating
**Post:** Titanium alloy
**Spherical Lock Nut and Adapter:** Cobalt Chrome alloy

### K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

### L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides will be present during routine MR imaging.

detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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