# English

## A. DEVICE DESCRIPTION

The Arthroereisis implant is a cannulated, tapered screw offered sterile in a range of sizes.

## B. INDICATIONS

### 1. Metal Implant Only

The Arthrex ProStop Arthroereisis Subtalar Implant is intended to assist in treating the hyperpronated foot by stabilizing the subtalar joint. It is intended to block forward, downward, and medial displacement of the talus thereby limiting excessive eversion of the hindfoot.

Examples include:

- Symptomatic acquired flat foot treatment in children, and adolescents
- Symptomatic congenital flexible flat foot
- Tarsal coalitions when associated with a flatfoot deformity
- Posterior tibial tendon dysfunction with supple feet
- Paralytic flat foot

### 2. Bioabsorbable Implant Only

The Arthrex ProStop Plus Arthroereisis Subtalar Implants are indicated as an internal support to primary surgical interventions in the treatment of flat foot, providing structural support at minimum during the first three months of healing.

## C. CONTRAINDICATIONS

1. Insufficient quantity or quality of bone.
2. Blood supply limitations and previous infections, which may retard healing.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
4. Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
5. Any active infection or blood supply limitations.
6. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
7. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
8. Do not use for surgeries other than those indicated.

## D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Foreign body reactions.
3. Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

## E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician
2. This device is intended to be used by a trained medical professional.
3. Factor in the relationship of the patient's weight to the size and material of the implant. Heavier patients require larger implants and may require a metal implant rather than a bioabsorbable implant.
4. An internal fixation device must never be re-used.
5. Do not re-sterilize this device.
6. All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
7. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.

8. Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
9. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
10. Detailed instructions on the use and limitations of the device should be given to the patient.
11. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
12. Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
13. Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
14. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
15. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

## F. MRI SAFETY INFORMATION

1. This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.

## G. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
2. Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

## H. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

## I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

## J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

This device is made of titanium or poly (L-lactide) (PLLA).

## K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-Bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

## L. INFORMATION

1. **In EU Only:** Procedures carried out using these devices may be used on the general population.
2. **In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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