A. DEVICE DESCRIPTION

Arthrex gravity tubing consists of two and four bag spike versions.

NOTE: All Arthrex gravity tubing is latex free.

B. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

C. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

D. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

E. PRECAUTION

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
2. Tubing systems manufactured in the EU do not contain phthalates, which are classified as Category 1 or Category 2 carcinogenic, genotoxic, and toxic substances under Annex I of Directive 67/548/EEC. Tubing systems manufactured in the USA may or may not contain phthalates. This ingredient is indicated on the packaging of the product.

F. DIRECTIONS FOR USE

WARNING: USING FLUID TO DISTEND ANY JOINT CARRIES WITH IT THE POS-

SIBILITY OF FLUID EXTRAVASATION INTO SURROUNDING TISSUE. Always use the lowest possible pressure settings to achieve the desired amount of distension and control of bleeding.

WARNING: PROPER OPERATION OF THIS DEVICE REQUIRES THE ESTABLISHMENT OF ADEQUATE FLUID OUTFLOW AND MONITORING OF THE SURGICAL FIELD. Carefully monitor and maintain fluid outflow and assess the patient regularly to avoid extravasation or any other adverse patient condition.

1. Prepare the arthroscope or the cannula for use with the tubing set.
2. Using sterile technique, carefully remove the tubing set from the package and pass it onto the sterile field.
3. Close all clamps, check the integrity of the remaining connections, and transfer the bag spikes out of the sterile field.
4. Using proper technique, insert the spikes into the fluid bags.

NOTE: When only one bag of liquid is to be connected, the spike on the unused bag must remain clamped.

5. Open all clamps.
6. Raise the fluid bags to the desired height.

NOTE: For every one foot of height above the patient, 22 mmHg of pressure is realized, or for every meter of height above the patient, 70 mmHg of pressure is realized.

7. To avoid air bubbles in the system, let fluid flow through the tubing until all air is purged.
8. Upon completion of the arthroscopic procedure, close all tubing clamps.
9. Dispose of the tubing in an acceptable waste container.

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