English

A. DEVICE DESCRIPTION

The Arthrex iBalance® TKA (Total Knee Arthroplasty) System consists of femoral, tibial tray, tibial bearing components and patellar components. All components are available in a range of sizes to fit varying anatomical requirements. Femoral components and tibial bearing components are available in both posterior stabilized (PS) and cruciate retaining (CR) configurations. Femoral components are available in left and right versions and are designed to work with the Arthrex dome patella component.

These devices are single use only and are intended for implantation with bone cement, with the exception of the porous coated femoral components which can be used cemented or uncemented (biological fixation).

B. INDICATIONS

The Arthrex iBalance TKA System is indicated for use in individuals undergoing surgery for:

- Painful, disabling joint disease of the knee resulting from degenerative arthritis, rheumatoid arthritis or post-traumatic arthritis
- Post-traumatic loss of knee joint configuration and function
- Moderate varus, valgus, or flexion deformity in which the ligamentous structures can be returned to adequate function and stability
- Revisions of previous unsuccessful knee replacement or other procedures

Additional indications for posteriorly stabilized components:

- Ligamentous instability requiring implant bearing surfaces with increased constraint
- Absent or non-functioning posterior cruciate ligament

These devices are single use only and are intended for implantation with bone cement, with the exception of the porous coated femoral components which can be used cemented or uncemented (biological fixation).

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- Do not use for surgeries other than those indicated.
- An overweight or obese patient can produce loads on the prosthesis which can lead to failure of fixation or failure of the device itself.
- Severe deformity and/or recurrent subluxation of the knee joint.
- Ligamentous instability or soft-tissue laxity such that the postoperative stability afforded by the prosthesis may be compromised or dislocate.

D. ADVERSE EFFECTS

- Loosening, migration, or loss of fixation of the implant.
- Foreign body reactions.
- Infections, both deep and superficial.
- Particulate wear debris and mid tissue discoloration from metallic components have been noted with other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
- Wear and damage to the articulating surfaces or soft tissue support structures.
- Postoperative symptoms include, but are not limited to: pain, persistent swelling, stiffness, limited range of motion, or incomplete resolution of preoperative symptoms.
- Periarticular calcification or ossification with or without impediment of joint mobility.
- Incomplete range of motion due to improper selection or positioning of components.
- Fatigue fracture of the implants.
- Intraoperative or postoperative bone fracture.
- Allergic reactions to materials.
- Other potential adverse effects of knee surgery may include: locking, neurovascular damage, dislocation, bending of the components, leg length discrepancies, and other less common adverse effects.
- Serious complications may be associated with any total joint replacement surgery. These complications include, but are not limited to genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Particles from metal and polyethylene components can be shed from the components during normal use and over time. These particles may be incorporated into the joint through natural biological processes or migrate to other parts of the body. The long-term effects of these particles may include, but are not limited to, cancer, lymphadenopathy and accumulation in other tissues and organs, and systemic disease.

E. WARNINGS

- Cautio: Federal law restricts this device to sale by or on the order of a physician.
- The Arthrex PS femoral components can function only with the Arthrex PS tibial bearing components. The Arthrex CR femoral components can only function with the Arthrex CR and CR Plus tibial bearing components. PS components should not be used with CR or CR Plus components.
- Care should be utilized in the handling of the components to minimize contamination and damage of the component surfaces.

- When using cement for fixation, the surgeon should use care to ensure complete cement contact and coverage on all parts of the prosthesis embedded in bone cement. Remove any excess cement to ensure that no cement remains on the articulating surfaces of the prosthesis.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- The surgeon must not allow damage to polished bearing surfaces because this may accelerate wear of the components. Any alteration or damage to a component may result in failure under load. Any prostheses so damaged must not be used.
- Components of the iBalance TKA System should not be used with those of another manufacturer since articular and dimensional compatibility cannot be assured.
- Detailed instructions on the use and limitations of this device should be given to the patient. Patient conditions can affect results and outcomes. Performing activities that increase stresses on implants such as running, lifting, skiing, etc. can result in early failure of these implants.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Postoperatively and until healing is complete, the fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- This product’s useful life will vary between patients, and is dependent upon accumulated weight-bearing activities and lifestyle.
- Increased likelihood of falls.
- Physicians should carefully assess the activity level of the patient prior to performing arthroplasty of the knee. Increased activity as well as increased weight can lead to accelerated wear of UHMWPE components.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex iBalance TKA is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only.
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system.
- Patient knee must be in direct contact with the patient table. Pillows should not be placed under the knee during scanning.

Under the scan conditions defined, the Arthrex iBalance TKA is expected to produce a maximum temperature rise of 2 °C at 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Arthrex iBalance TKA extends approximately 20 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Proper selection of the Arthrex iBalance TKA System depends on the judgment of the surgeon.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- The femoral, tibial tray and tibial bearing components are supplied individually in a double layer packaging configuration. The components should not be used if the seal or packaging has been or is suspected of having been breached.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Femoral and Tibial tray: Cobalt Chrome and/or Titanium
Tibial bearing and patellar components: Ultra high molecular weight polyethylene (UHMWPE)

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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