A. DEVICE DESCRIPTION

The Arthrex AnaToemic® Phalangeal Prosthesis is a one piece implant system. The implant has a polished concave oval disk and a stem with barbs for fixation in bone. The prosthesis is available in a range of sizes to fit varying anatomical requirements.

B. INDICATIONS

The Arthrex AnaToemic Phalangeal Prosthesis is a press-fit implant that is intended to be used in patients with hallux rigidus, hallux limitus, hallux valgus, arthritic degradation of the metatarsophalangeal joint, degenerative arthritis, rheumatoid arthritis, and bunion deformity associated with arthritis of the metatarsal-phalangeal joint.

C. CONTRAINDICATIONS

1. Insufficient quantity or quality of bone.
2. Severe tendon, neurological, or vascular deficiencies related to the lower extremities.
3. Blood supply limitations and previous infections, which may retard healing.
4. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
5. Any active infection or blood supply limitations.
6. The prosthesis should not be used in cases of severe Hallux Valgus deformity or in combination with a metatarsal or phalangeal osteotomy.
7. Do not use for surgeries other than those indicated.
8. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
9. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Foreign-body reactions.
3. Limited motion up to stiffness of the joint.
4. Erosion of the cartilage of the metatarsal head.
5. Postoperative Hallux deformity.
6. Resistant pain.

E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician
2. This device is intended to be used by a trained medical professional.
3. Do not re-sterilize this device.
4. All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
5. An internal fixation device must never be reused.
6. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
7. Detailed instructions on the use and limitations of the device, the patient leaflet and the patient implant card should be given to the patient.
8. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should

be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.

9. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.

10. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.

11. **TPS/CaP and TiNBN coated device - Not available in the US:** Fluid contact other than patient's blood should be avoided to achieve best ongrowth results.

12. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.

13. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION**1. MR Conditional**

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the AnaToemic Phalangeal Prosthesis are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the AnaToemic Phalangeal Prosthesis is expected to produce a maximum temperature rise of 4.6 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the AnaToemic Phalangeal Prosthesis extends approximately 16 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
2. Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

H. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for material specifications.

This device is made of Cobalt Chrome. AnaToemic prostheses are available uncoated and with TiNBN or TPS/CaP coatings.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

1. **In EU Only:** Procedures carried out using these devices may be used on the general population.
2. **In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

M. POSTOPERATIVE REHABILITATION PROTOCOL

The progression to normal ambulation and the use of standard foot gear is limited only by the persistence of postoperative swelling and discomfort.

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