A. DEVICE DESCRIPTION

The Arthrex iBalance® High Tibial Opening Wedge Osteotomy (HTO) System is comprised of an anatomically designed implant, fixation anchors, and a surgical instrument system. The wedge-like implants, and cortical and cancellous fixation anchors are available in various sizes.

B. INDICATIONS

The Arthrex iBalance High Tibial Opening Wedge Osteotomy (HTO) System includes instrumentation, implant, and fixation and is used to provide fixation following proximal tibial opening wedge osteotomies. It is intended to be used with adequate post-operative immobilization.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.
- Patients who have not reached skeletal maturity.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be reused. The devices may be contaminated and not properly sterilized. This may result in infection, toxic reaction or similar adverse effect. Re-use of a device may result in mechanical damage such as deformation of the threads on the implant or anchors, weakening of the anchor, or deformation of the feature that inhibits back-out of the anchor. This may cause difficulties in implantation or failure of the device post-operatively.
- Do not re-sterilize this device.
- Post-operatively, the patient is placed in a post-op, hinged knee brace locked in 0 degrees of extension.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex iBalance instrument system is required for proper implantation of the device. Weight bearing should not begin until there is evidence of bone growth.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate post-operative management.

- Detailed instructions on the use and limitations of the device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the iBalance High Tibial Opening Wedge Osteotomy System is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 1-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the iBalance High Tibial Opening Wedge Osteotomy System is expected to produce a maximum temperature rise of up to 6 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the iBalance High Tibial Opening Wedge Osteotomy System can extend up to approximately 120 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Before opening the implant sterile barrier pouch, confirm implant size by viewing the engraving on the implant through the clear packaging.
 - Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
 - Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- CAUTION:** Improper alignment of the instrumentation may result in unintended tibial A/P slope change.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The devices are made out of polyetheretherketone (PEEK) and tantalum.

K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-Bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

- In EU Only:** Procedures carried out using these devices may be used on the general population.
- In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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