English

A. DESCRIPTION

The Apollo^{RF} probes are designed specifically for use with the Synergy^{RF}™ (Radio Frequency) Console only and are not compatible with any other electrological generator.

Non-Aspirating RF probes provide a bipolar electrological effect to target tissue. Located on the top of the RF probe handle, the black button adjusts the default power setting for the specific RF probe connected to the Synergy^{RF} Console. The yellow button invokes the Ablate function. The blue button that is farthest from the end of the RF probe invokes the Coagulation function. A Synergy^{RF} Footswitch is connected to the front panel of the Synergy^{RF} Console to override the RF probe’s control buttons. The user has the option to turn off the footswitch override through the Synergy^{RF} Console touch screen options. **Aspirating RF probes** provide a bipolar electrological effect with control buttons as described above and can aspirate materials away from the surgical field.

B. INDICATIONS

The RF probe are accessories to the Synergy^{RF} Console and are intended for use as a complete system in the resection, ablation, and coagulation of soft tissue and hemostasis of blood vessels and tissue in arthroscopic and orthopedic procedures. Specifically, the RF probes, Synergy^{RF} Console and their accessories are used for arthroscopic surgery of the shoulder, wrist, hand, elbow, hip, knee, foot and ankle.

C. CONTRAINDICATIONS

- Use of an RF probe in any non-arthroscopic surgical procedure.
- Use of an RF probe during an arthroscopic procedure without a conductive irrigation solution.
- Use of an RF probe in a surgical procedure for a patient where an arthroscopic procedure is contraindicated for any reason.
- Use of an RF probe in patients with cardiac pacemakers or other electronic device implants.

D. ADVERSE EFFECTS

- As a consequence of electrography, damage to surrounding tissue through iatrogenic injury could occur.
- Additional damage to tissue from excessive force during blunt dissection is possible.
- Premature tip wear may be caused by vigorous use against bony surfaces.
- Thermal damage to joint tissue or dermal burns around portal incisions are possible, if fluid flow is not continuous.

E. PRECAUTIONS

- Prior to initial use, ensure that all package inserts, including warnings, precautions, *User’s Guides* and Directions for Use are read and understood.
- The RF probes are supplied sterile and intended for SINGLE USE ONLY. DO NOT clean, sterilize, or reuse the RF probe.
- DO NOT use excessive force when inserting or removing the RF probe. DO NOT insert the RF probe into an obstructed passageway. Patient injury and or product damage may occur.
- DO NOT use an RF probe as a lever to enlarge the surgical site or gain access to tissue, which may result in a bent or damaged RF probe.
- Always keep the tip of the active RF probe completely immersed in conductive irrigation fluid (saline, Ringers’ lactate, etc.), regardless of the tpe of arthroscopic surgery being performed. DO NOT use pre-warmed conductive irrigation fluids as this may result in tissue damage or thermal injury.
- A continuous flow of irrigation fluid is necessary during use of the Apollo probes. Fluid flow assists in removing tissue debris and reducing the intra-joint temperature between activations. RF probes utilized in stagnant fluid can result in heated irrigant in the joint, which can result in tissue damage inside the joint, skin burns near portals or result in damage to the probe.
- For RF probes, failure to attach the Suction Adapter to an external suction source may cause thermal injury to the patient or user.
- The rate and depth of tissue ablation can be affected by a number of contributors, such as, but not limited to the style of RF probe selected, the Synergy^{RF} Console power setting, amount of pressure exerted on the target tissue by the RF probe tip, and the speed that the RF probe is moved over the target tissue to be ablated.
- DO NOT place cables from other equipment between the probe and Synergy^{RF} Console.
- DO NOT allow patient contact with grounded metal objects.
- Prolonged ablation should be avoided. Intermittent pauses in ablation are recommended, to allow warm fluid and tissue debris to be evacuated from the joint. Prolonged probe activation while using the suction feature may result in elevated shaft or suction tubing temperatures.

F. WARNINGS

- Be sure to read and understand the Synergy^{RF} Console *User’s Guide (Arthrex DFU-0221-XX)* prior to using any RF probes as they are designed to be operated as a complete system.
- DO NOT use an RF probes in the presence of flammable materials. Do not place activated RF probes near or in contact with flammable materials (gauze, surgical drapes, etc.). Sparking and heating associated with electrography may be an ignition source and fire could result.
- DO NOT use non-conductive irrigation fluids (sterile water, air, gas, glycine, Purisole, sorbitol-mannitol, etc.).
- DO NOT touch the tip of the RF probe while activated.
- DO NOT contact metal objects while the RF probe is activated, as this may cause damage to the RF probe.

probe slightly off of the tissue to be treated. It is always good practice to keep the tip of the RF probe moving along, slightly off the tissue being treated.

- WARNING:** An Active RF probe tip must never be outside of the user’s field-of-view during surgery, as this can result in unintentional damage to tissue.
- NOTE:** It is normal to view bubbles around the shaft of the RF probe as well as an illumination of its tip when active.
- NOTE:** When a Synergy^{RF} Footswitch is connected to the front panel of a Synergy^{RF} Console, its controls override the buttons of the RF probe’s handle. The user has the option of overriding this feature using the Synergy^{RF} Console touch screen. Please refer to the Synergy^{RF} Console *User’s Guide (Arthrex DFU-0221-XX)* for complete information.
- To use the Ablate function:** Press the yellow Ablate button on the top of the RF probe’s handle to invoke the Ablate function, or if using a Synergy^{RF} Footswitch, depress the yellow pedal on the left side of the Synergy^{RF} Footswitch.
- To use the Coagulation function:** Press the blue Coag button on the top of the RF probe’s handle to invoke the Coagulation function, or if using an Synergy^{RF} Footswitch, depress the blue pedal on the right side of the Synergy^{RF} Footswitch.
- To use aspirating RF probes:** After connecting to a suction source, open the roller clamp of an aspirating RF probe fully prior to and during use. To shut off the aspiration feature, close the clamp on the suction tube only when the RF probe is outside of the patient.
- To remove eschar from the RF probe tip:** Remove the RF probe from the patient, rub the RF probe tip against an abrasive pad, and then against a wet gauze pad.

I.V. TERMINATING SURGERY

- Push the main power switch to the “OFF” position on the front panel of the Synergy^{RF} Console to deactivate the Synergy^{RF} System.
- Disconnect the RF probe’s connector from the front panel of the Synergy^{RF} Console.
- Disconnect the suction tubing.
- Dispose of the RF probe by following standard hospital protocol.

L. ALARMS AND TROUBLESHOOTING

The Synergy^{RF} Console issues an audible continuous tone whenever the RF probe is active. It also issues an audible alarm when it detects a significant change in resistance. Refer to the Synergy^{RF} Console *User’s Guide (Arthrex DFU-0221-XX)* or contact your Arthrex representative for more information.

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