English

A. DEVICE DESCRIPTION

The Arthrex Corkscrew™, PushLock®, and SwiveLock® suture anchors consist of anchors with an integral or separate eyelet. The PushLock™ Tenodesis anchor is a two-piece insert-in anchor with either a forked or closed eyelet. They are pre-loaded on a handled inserter. Suture, with or without needles, and a suture threader may also be provided.

The Arthrex Implant System, Tenodesis Screw Eyelet is a kit comprised of attachable eyelets, FiberTape™ suture and a suture threading device. This Implant System is used in conjunction with Arthrex Tenodesis Screws and a Tenodesis driver for use as a SwiveLock (Tenodesis) suture anchor.

B. INDICATIONS

The Arthrex Corkscrew, PushLock, and SwiveLock suture anchors are intended to be used for suture (soft tissue) fixation to bone in the foot, ankle, knee, hand, wrist, elbow, shoulder, and hip (devices with FiberWire® only): Acetabular labral repair (**except** DX SwiveLock SL, 2.4 mm and 2.5 mm PushLock, Nano, Micro and Mini Corkscrew FT suture anchors). The Arthrex PushLock Tenodesis anchor is intended to provide soft tissue re-attachment (i.e. fixation of ligament and tendon graft tissue) in surgeries of the shoulder, elbow, knee, foot/ankle, and hand/wrist.

Surgeons must apply their professional judgment when determining the appropriate suture anchor size based on the specific indication, preferred surgical technique, and patient history:

Shoulder: Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Debrid Repair, Capsular Shift or Capsulolabral Reconstruction

Foot/Ankle: Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Mid-foot Reconstruction, Metatarsal Ligament Repair, Metatarsal Tendon Repair (**except** PushLock Tenodesis anchor) and Bunionectomy (**except** PushLock Tenodesis anchor, 2.4 mm and 2.5 mm PushLocks, DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors). **DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors only:** Digital tendon transfers. **PushLock Tenodesis anchor only:** Flexor Hallucis Longus for Achilles Tendon reconstruction, tendon transfers in the foot and ankle

Knee: Anterior Cruciate Ligament Repair (**except** PushLock, SwiveLock, DX SwiveLock, DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors), Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair and Iliotibial Band Tenodesis

Hand/Wrist: Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction. **2.5 mm PushLock, DX SwiveLock SL, Nano, Micro and Mini Corkscrew suture anchors only:** Repair/Reconstruction of Collateral Ligaments, Carpal Ligament Reconstruction, Repair of Flexor and Extensor Tendons at the PIP DIPr and MCP joints for all digits, and Digital Tendon Transfers. **PushLock Tenodesis anchor, DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors only:** Carponmetacarpal Joint Arthroplasty (Basal Thumb joint Arthroplasty). **PushLock Tenodesis anchor only:** Carpal Ligament Reconstructions and repairs, tendon transfer in the hand/wrist

Elbow: Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction, Lateral Epicondylitis repair (**except** PushLock Tenodesis anchor, 2.4 mm and 2.5 mm PushLocks, DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors)

The 2.4 mm Hip PushLock suture anchor is indicated for acetabular labral repair **ONLY**.

Gluteus Medius (U.S. Only): 4.75 – 5.5 mm PEEK and Titanium Corkscrew suture anchors and 5.5- 6.5 mm PEEK and Titanium Corkscrew suture anchors

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Bioabsorbable only:** Foreign Body Reactions. See Adverse Effects- Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Bioabsorbable only:** Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Titanium only:** Shoulder dislocation/subluxation.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be reused.
- Do not re-sterilize this device.
- Titanium only:** All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Bioabsorbable only:** Attempting implantation into hard, dense bone and/or drilling/punching smaller diameter holes than recommended may cause failure (breakage) of the implant during insertion.
- Bioabsorbable Corkscrew suture anchor only:** The Arthrex 6.5 mm anchor should be used in soft bone only.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this

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