A. DEVICE DESCRIPTION

The Arthrex Autologous Conditioned Plasma (ACP®) Double Syringe System (Arthrex ACP System) is comprised of an inner syringe inside an outer syringe. The inner syringe is specifically designed to fit within and function in concert with the specifically designed outer syringe. The two syringes are connected via a luer lock connector. A threaded screw cap connected to the luer tip of the outer syringe can be removed prior to use. The Arthrex ACP System was designed to be used with benchtop centrifuges with swing-out rotors that securely accommodate ACP Double Syringes (e.g., Roto-fix 32 A Benchtop Centrifuge).

B. INDICATIONS

For the USA, the devices comprised of the Arthrex ACP System are indicated for:

- The Double Syringe Autologous Conditioned Plasma (ACP) System is used to facilitate the safe and rapid preparation of autologous platelet-rich plasma (PRP) from a small sample of blood at the patient's point of care.

The PRP can be mixed with autograft and allograft bone prior to application to an orthopedic surgical site as deemed necessary by the clinical use requirements.

For outside the USA the devices comprised of the Arthrex ACP System are indicated for:

The Autologous Conditioned Plasma (ACP) Double Syringe System is used for sterile separation of non-homogenous liquids. This device:

- Will allow the withdraw of fluids from the human body using a commercially available cannula with LuerLock connection;
- May be centrifuged up to RCF 1100g thereafter to transfer the supernatant into the smaller syringe by pressure balance.

C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Allergies and other reactions to device materials.
3. Hematoma.
4. Damage to blood vessels and nerve damage resulting in pain or numbness from autologous sampling.
5. Delayed wound healing.

E. WARNINGS

1. Do not re-sterilize this device.
2. Follow the manufacturer's instructions when using centrifuge. Use only benchtop centrifuges with swing-out rotors that securely accommodate ACP Double Syringes (e.g., Roto-fix 32 A Benchtop Centrifuge). Outcomes using other types of centrifuges are unknown.
3. Follow the manufacturer's directions for use if using an anticoagulant (e.g., Anticoagulant Citrate Dextrose Solution A (ACD-A)).
4. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised me-

chanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

5. Failure to use this device in accordance with the directions for use may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.

F. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

H. STERILIZATION

The device is supplied sterile. Refer to the package label for the sterilization method.

I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

- Syringe and Piston:** Polypropylene, TPE (Thermoplastic Elastomer), silicone fluid
- Cap:** Polyethylene HD (High Density)
- Small Syringe:** Polypropylene, Polyisoprene rubber, silicone fluid

J. STORAGE CONDITIONS

These devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

K. INFORMATION

Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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