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AR-6572S

2.6 mm diameter  
x 305 mm length

for use with the  
Gemini 8.25 mm  
Self-Retaining  
Cannula System

Manufactured for:

**Arthrex, Inc.**

1370 Creekside Blvd.

Naples, FL 34108 • USA

1-(800) 934-4404

www.arthrex.com

EC | REP

**Arthrex GmbH**

Erwin-Hielscher-Strasse 9

81249 München, Germany

Tel: +49 89 909005-0

www.arthrex.de

**Frantz Medical Development**

7740 Metric Drive

Mentor, OH 44060, USA

(440) 255-4514

©2013 Frantz Medical Ltd.  
All rights reserved.

English	201004B 4/15
---------	--------------

Please read all information carefully. Failure to follow the instructions may lead to user and/or patient injury.

NOTE: This insert is intended to provide instructions for the Switching Stick for use with Gemini Self-Retaining Cannula System. It is not a reference for surgical techniques.

#### A. SUPPLIED

The Gemini Switching Stick is provided non-sterile for multiple use.

#### B. INDICATIONS

The Switching Stick is a reusable instrument used as a guide with the Gemini 8.25 mm Self-Retaining Cannula System.

#### C. CONTRAINDICATIONS

The Switching Stick is intended for use only as indicated.

#### D. CLEANING AND STERILIZATION PROCEDURE

The following procedure is a recommendation to be used in conjunction with currently accepted professional standards for cleaning and sterilization of the Switching Stick. Frantz Medical assumes no responsibility for providing parameters for the cleaning and sterilization of the device, for user negligence associated with the cleaning and sterilization procedure, or for intentional misuse of the device.

#### E. CLEANING

It is recommended that the switching stick undergo cleaning as soon as possible after use; immediate removal of excess soil with disposable wipes will prevent drying of adherent bodily fluids and debris.

- Manual Cleaning:**

Scrub the Switching Stick with a soft bristle brush in warm water (20 - 25°C/68 - 77°F) with a mild detergent/enzymatic cleaner. Thoroughly rinse the device with clean water after the detergent/enzymatic wash.
- Automatic Cleaning / Disinfection:**

The Switching Stick may be placed in an automatic dishwasher at 90°C (194°F).

#### F. STERILIZATION

The device must be sterilized prior to surgical use. Sterilization procedures should only be conducted by persons with proper knowledge and training techniques.

Steam sterilization is recommended using the Association for the Advancement of Medical Instrumentation (AAMI) guideline “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.” The cycle parameters recommended in this standard are suggestions. Cycle parameters should always be verified against the sterilizer manufacturer’s instructions for the specific sterilizer and load configuration being used.

The device must be adequately cooled, after being removed from the sterilizer. Do not touch the device during cooling or attempt to facilitate the cooling process by placing the device on a cold surface or immersing it in cold fluid.

#### G. WARNINGS AND PRECAUTIONS

- Arthroscopic procedures should only be conducted by persons with proper knowledge and training techniques. Prior to a minimally invasive, arthroscopic procedure, it is the responsibility of the person performing the procedure to be familiar with all medical literature, methods, and potential complications and hazards that may be associated with such a procedure.
- The Switching Stick is provided **non-sterile**. All Switching Sticks must be cleaned and sterilized prior to initial use as well as reuse.
- When properly performed, cleaning and sterilization of the Switching Stick will not compromise the performance or integrity of the device. Inspect the device for damage at all stages of handling. *If physical damage is observed, consult the manufacturer for guidance.*
- Detergents and solutions should have a pH between 7.0 and 9.5.
- Sterilization temperature should not exceed 135°C (275°F).

Deutsch	201004B 4/15
---------	--------------

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