A. DEVICE DESCRIPTION

The RetroDrill® Guide Assembly is a drill with graduated markings. It's available as cannulated, or non-cannulated.

RetroCutters® have reverse flutes to cut from the inside – out, and are available in many sizes. Dual RetroCutters have cutting surfaces on both sides, and are available in many sizes.

B. INDICATIONS

The Arthrex RetroDrill Guide Assembly is intended to assist in the creation of bone sockets and tunnels. The reverse-cutting flutes of the Arthrex RetroCutter are intended to allow retrograde socket and tunnel creation with the use of the RetroDrill Guide Pin. Dual RetroCutters have both reverse and forward cutting flutes for retrograde and antegrade socket and tunnel creation. The RetroDrill Guide Pin is circumferentially laser-marked in 5mm increments to assist in measuring tunnel depth.

C. CONTRAINDICATIONS

1. Insufficient quantity or quality of bone.
2. Blood supply limitations and previous infections, which may retard healing.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to surgery.
4. Any active infection or blood supply limitations.
5. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
6. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
7. Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Foreign body reactions.

E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician
2. This device is intended to be used by a trained medical professional.
3. Do not re-sterilize this device.
4. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
5. To prepare to drill the socket in bone, first verify that the directional switch on the drill is set to the FORWARD (clockwise) direction. Then, simultaneously activate the drill while pulling back to create the bone socket.
6. When reloading the RetroCutter back onto the RetroDrill Guide, assure that the RetroDrill Guide Pin tip is engaged into the RetroDrill Guide before drilling in reverse (counterclockwise).
7. Always match the Adapteur™ Guide C-Ring angle with the angle of the chosen RetroDrill Guide and manually test the alignment with the RetroDrill Guide Pin before use with a power drill.
8. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to

wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

9. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
10. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
2. Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

H. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The RetroDrillGuide Pin and RetroCutter are made of stainless steel.

The O-Ring is made of Viton® fluorocarbon elastomer.

The Nitinol Guide Pin is made of nitinol.

J. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

K. INFORMATION

1. **In EU Only:** Procedures carried out using these devices may be used on the general population.
2. **In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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