Symbols Used on Labeling
Auf der Verpackung verwendete Symbole
Symboles utilisés sur l'étiquetage
Simboli usati sull'etichetta
Símbolos utilizados en las etiquetas
Símbolos utilizados na rotulagem

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English **DFU-0243-2r0**

A. DEVICE DESCRIPTION
The BioSync® Reconstruction Wedge is a porous metal wedge used for angular correction of small bones in the ankle and foot. It is offered with varying widths and thicknesses to accommodate a variety of small bone applications.

B. INDICATIONS
The BioSync Reconstruction Wedge is intended to be used for internal bone fixation for bone fractures, fusions, or osteotomies in the ankle and foot, such as:

- Cotton and Evans Wedges:**
1. Opening wedge osteotomies of the bones of the foot (including osteotomies for Hallux Valgus, except in Canada)
 2. Opening wedge of Medial Cuneiform or Cotton osteotomies
 3. Lateral Colum Lengthening (Evans Lengthening Osteotomy or Calcaneal Z Osteotomy)
 4. Metatarsal/Cuneiform arthrodesis

Midfoot Wedges:

1. Opening wedge osteotomies of the bones of the foot including osteotomies for Hallux Valgus
2. Nonunion of arthrodesis of the Midfoot including Metatarsal/Cuneiform arthrodesis (TMT or Lapidus)

This device is intended for use with ancillary fixation. The BioSync Reconstruction Wedge is not intended for use in the spine.

- C. CONTRAINDICATIONS**
1. Insufficient quantity or quality of bone.
 2. Blood supply limitations and previous infections, which may retard healing.
 3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
 4. Any active infection or blood supply limitations.
 5. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
 6. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
 7. Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Foreign body reactions.
3. Unwanted soft tissue damage and / or articular surface damage due to improper wedge placement.

- E. WARNINGS**
1. An internal fixation device must never be re-used.
 2. Do not re-sterilize this device. All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
 3. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
 4. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
 5. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
 6. Detailed instructions on the use and limitations of this device should be given to the patient.
 7. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
 8. Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

9. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.

F. MRI SAFETY INFORMATION

1. This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.

G. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.

I. STERILIZATION
This device is provided sterile. Refer to the package label for the sterilization method.
Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information. .

J. MATERIAL SPECIFICATIONS
Refer to the package label for the materials.
This device is made of titanium.

K. STORAGE CONDITIONS
This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION
Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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