<p><b>Symbols Used on Labeling</b></p> <p>Auf der Verpackung verwendete Symbole</p> <p>Symboles utilisés sur l'étiquetage</p> <p>Simboli usati sull'etichetta</p> <p>Símbolos utilizados en las etiquetas</p> <p>Símbolos utilizados na embalagem</p>
---

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## English DFU-0203r I

### A. DESCRIPTION

Arthrex Footswitches are designed specifically to work with Arthrex consoles, be sure to match the correct footswitch to the correct console. AR-6483 and AR-6483W are designed to control Arthrex DualWave Arthroscopy Pump (AR-6480). AR-8310 and AR-8310W are designed to control Arthrex Adapteur Power System™ II (APS II) (AR-8300) or Arthrex Synergy™ Resection Console (AR-8305). AR-8310SP is designed to control both the Arthrex DualWave Arthroscopy Pump (AR-6480) and the Arthrex Adapteur Power System™ II (APS II) (AR-8300) or Arthrex Synergy™ Resection Console (AR-8305).

The footswitch lets you adjust the settings at a distance from the control unit of the console.

For more details on the effect of each function refer to console's manual that the footswitch is being connected to.

Note: Each footswitch is designed to control specific functions for each console.

### B. CAUTIONS

- Do NOT sterilize this device.
- Do NOT autoclave this device.
- Do NOT clean and disinfect in a thermo washer disinfector.
- NEVER immerse this device in water or cleaning detergent.
- Do NOT clean the device with abrasive cleaning or disinfectant compounds, solvents, or other materials that could scratch or damage the device.
- Do NOT lift the device by the cord.
- Keep this device away from magnetic and electromagnetic fields.
- This unit is NOT suitable for use in the presence of flammable mixtures.

### C. CLEANING AND DISINFECTION

- The footswitch is supplied non-sterile and is not to be sterilized.
- To clean and disinfect, unplug the cord from the console and add the protective cap to the shaver connector.
- Use a commercially available cleaning disinfectant for surfaces.
- Alkaline agents with a pH value > 8 are to be avoided.
- Always comply with the instructions issued by the manufacturer of the cleaning disinfectant regarding concentration, exposure times, temperature and material compatibility.
- Do NOT clean the device with abrasive cleaning or disinfectant compounds, solvents, or other materials that could scratch or damage the device.
- The footswitch is NOT suitable to be cleaned and disinfected in a thermo washer disinfector.  
CAUTION: Do NOT allow fluids to enter the connector or allow the receptacles to have any contact with liquids. Remove dust or moisture, if necessary, with dry compressed air. ONLY dry connectors should be plugged into the console.

### D. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact customer service if the package has been opened or altered.

### E. INSPECTION

Before use, the device should be inspected for signs of wear or damage. In particular look for any of the following problems:

- Damage to the cord
- Damage to the cord receptacle
- Damage to the connector
- Damage to the buttons
- Obvious damage to the unit
- Traces of blood, debris, and stains

For continued safety and performance of the footswitch, it is recommended that the functionality and electrical safety of the switch to be verified on an annual basis by a qualified person.

### F. STORAGE CONDITIONS

The device should be stored at room temperature in a dry environment.

### G. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical techniques information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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