English

A. DEVICE DESCRIPTION

The Biceps Tenodesis Implant Systems are available in multiple configurations for use in biceps repair. One configuration consists of a BicepsButton™, a Tenodesis Screw (PEEK or BioComposite), suture and manual surgical instruments. Another configuration consists of a Tenodesis Button, suture and manual surgical instruments. Another configuration consists of Pec button assemblies and manual surgical instruments. Another configuration consists of a knotted button assembly, suture, and manual surgical instruments. Another configuration consists of a SwiveLock® suture anchor (PEEK), suture and manual surgical instruments.

B. INDICATIONS

The BicepsButton, Pec Button and Tenodesis Button are used for fixation of bone to bone or soft tissue to bone, and is intended as fixation posts, a distribution bridge, or for distributing suture tension over areas of ligament or tendon repair in the shoulder and elbow. Arthrex recommends the use of #2 FiberWire® suture or equivalent with the BicepsButton, Pec Button and Tenodesis Button. The Arthrex Suture Family is intended for use in soft tissue approximation and or ligation. These sutures may be incorporated, as components, into surgeries where constructs including those with allograft or autograft tissues are used for repair.

The Tenodesis Screw is intended to provide soft tissue reattachment (i.e. fixation of ligament and tendon graft tissue).

The Arthrex SwiveLock suture anchors are intended to be used for suture (soft tissue) fixation to bone in the shoulder.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.
- Arthrex suture is not for use in cardiac indications.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Allergies and other reactions to device materials.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Silicone sensitivity, though very rare, has been reported.
- Adverse reactions have not been noted with the Arthrex suture products in animal testing.
- Common non-absorbable suture reactions may include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site.
- If needles are included, inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- Metal implants only:** All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device, the patient leaflet and the patient implant card should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Once open, discard unused suture.
- Do not expose suture to heat.
- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing

Arthrex FiberWire or Arthrex UHMWPE suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

- As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formations. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.
- Metal Implants Only:** Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional Information for Suture Buttons

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Suture Buttons manufactured from Titanium are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Suture Buttons are expected to produce a maximum temperature rise of up to 1.7° C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Suture Buttons can extend up to approximately 17 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

2. MR Safe Information

Arthrex Suture manufactured from UHMWPE and polyester with or without silicone elastomer coating, cyanoacrylate, and nylon are MR safe.

The Tenodesis Screws devices manufactured from only polyetheretherketone (PEEK), polypropylene or poly (L-Lactide), poly (L, DL-lactide; PLDLA), biphasic calcium phsphate, biphasic tricalcium phosphate (TCP) and/or hydroxyapatite (HA) are MR safe.

The SwiveLock suture anchor devices manufactured from only polyetheretherketone (PEEK), poly (L-Lactide, PLLA), Poly (L-Lactide-co-D acid, PLDLA), and/or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP) are MR safe.

G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate screw size based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Loop and button should be properly oriented during passing and fixation, per technique guide, in order to function properly. In handling suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- Assure that all knots have been secured using accepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Care should be taken to prevent damage to surrounding tissue or user puncture due to improper handling of the needlepoint.
- Do not grasp the needle if included at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.
- Under insertion of the device may leave the proximal end of the implant protruding beyond the cortical bone, which could potentially cause soft tissue irritation and/or pain post-operatively.
- Bio-Tenodesis™ Screw only:** Use the appropriate size Arthrex Drill to create a pilot hole in the bone.
- Make sure to use the recommended drill bit or punch to create the bone socket.
- SwiveLock suture anchors only:** During anchor insertion, ensure that the angle of anchor insertion is coaxial to that of the previously prepared bone socket.
- SwiveLock suture anchors only:** Insert the driver into the bone socket until the anchor body makes contact with the bone. Preview and adjust suture tension, if necessary. Tension will not increase during final advancement of the anchor body.
- SwiveLock suture anchors only:** Ensure that the anchor body is in full contact with the bone before advancing the anchor body into the prepared bone socket.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST179 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Button: This device is made of titanium.

Suture: Ultra High Molecular Weight Polyethylene (UHMWPE) braid, or a polyblend of UHMWPE and Polyester. The sutures meet or exceed U.S. and European Pharmacopaeia standards for non-absorbable surgical sutures (except for diameter requirements).

Additional materials to the sutures may include silicone elastomer coating, cyanoacrylate and/or nylon. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue. The suture is available non-dyed, dyed, and fully or partially striped. Suture dyes include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon.

Arthrex sutures, when tested per ISO 10993, Biological Evaluation of Medical Devices-Part 10: Tests for Irritations and Sensitization, had no reactions of allergic or sensitive nature. The dyed suture and coating are pharmacologically inactive.

Arthrex FiberWire suture is not absorbed, but may become encapsulated in the surrounding connective tissues. The Arthrex FiberWire is not known to have significant change in tensile strength in vivo.

Tenodesis Screws: polyetheretherketone (PEEK), polypropylene or poly (L-Lactide), poly (L, DL-lactide; PLDLA), biphasic calcium phosphate, biphasic tricalcium phosphate (TCP) and/or hydroxyapatite (HA)

SwiveLock suture anchors consist of either one or two components. Each component is manufactured from either/or polyetheretherketone (PEEK), Poly (L-Lactide acid, PLLA), Poly (L-Lactide-co-D acid, PLDLA), or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP).

K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C/ 90° F, and should not be used after the expiration date.

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

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