Symbols Used on Labeling
Auf der Verpackung verwendete Symbole
Symboles utilisés sur l'étiquetage
Simboli usati sull'etichetta
Símbolos utilizados en las etiquetas
Símbolos utilizados na rotulagem

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English

A. DEVICE DESCRIPTION
The PEEKPower™ Plates are available in different shapes, sizes and orientations (e.g. left and right types). The plates have specific sized holes for screws to provide fixation.

B. INDICATIONS
The PEEKPower Plates are used as an internal fixation device to stabilize fractures or osteotomies at the time of healing. The appropriate plate should be chosen for the site to be treated.

The PEEKPower Distal Radius Plate is intended for fixation of fractures of the radius.

The PEEKPower Humeral Fracture Plate is intended for fixation of proximal fractures of the humerus.

The PEEKPower Tibial Head Fracture Plates are intended for fixation of proximal fractures of the tibia.

The PEEKPower Osteotomy Plates provide fixation following opening wedge osteotomies.

The plates are intended to be used with adequate post-operative immobilization.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.
- PEEKPower High Tibial Osteotomy (HTO) and Lateral Distal Femoral Osteotomy (LDFO) Plates:** Previous closing wedge osteotomy in the same region of the same extremity.
- PEEKPower HTO and LDFO Plates:** Closing Wedge Osteotomy in the same region of the same extremity.
- PEEKPower HTO and LDFO Plates:** Rotational Osteotomy in the same region of the same extremity.

D. ADVERSE EFFECTS

- Foreign body reactions
- Infections, both deep and superficial.

E. WARNINGS

- Use only the Arthrex Screws provided with the PEEKPower Plate systems.
- Osteotomy plates only:** Corresponding preoperative planning must be followed.
- Postoperatively and until bone healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Removal of the device should be performed using standard surgical practices for device removal. Implant removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Do not resterilize this device.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- The PEEKPower Plates should not be used in patients whose body weight is more than 100 kg. (220 lbs.).
- An internal fixation device must never be reused.

F. MRI SAFETY INFORMATION

- This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment.

This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.

G. PRECAUTIONS

- The plate cannot be bent or intra-operatively contoured. Contouring the plate could result in premature failure/breakage.
- The locking screw fixation will be compromised if the plate is damaged during the drilling process.
- The locking screws can be inserted into the same plate hole a second time, at the same or different angle, if necessary. For example, if a locking screw needs to be angled differently in the hole, or a different length screw needs to be used.
- The screws should be inserted by hand with the appropriate screwdriver and not with powered instruments. The most effective fixation is achieved when the screw head is inserted flush with the top surface of the plate.
- Proper placement of the plate is critical to avoid impingement and to ensure optimal screw placement.
- PEEKPower, HTO and LDFO plate: For a fracture of the far cortex, a compression screw must be used to reduce the fracture fragments. Once the fragments are stabilized and reduced, the compression screws must be replaced with locking screws for the most stable fixation.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.

I. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.
PEEKPower Plates are made of Carbon Fiber Reinforced PEEK (CF-PEEK) and tantalum.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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PEEKPower™ Plates

PEEKPower Platten
Placas PEEKPower
Placas PEEKPower
Plaques PEEKPower
Piastre PEEKPower

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Tel: +49 89 909005-0

www.arthrex.de

This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the "Arthrex U.S. Product Warranty" section of the Arthrex, Inc. website, found at www.arthrex.com whose provisions are incorporated herein by reference.

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Arthrex GmbH