English

A. DEVICE DESCRIPTION

The PIP (Proximal Interphalangeal) Dart is an implant with a barbed design and ranges in size from 20 to 35 mm in length. The PIP Dart is inserted between the proximal and middle phalanges, so the barbs fixate on the phalangeal canal of the toe.

The RetroFusion™ screws is an implant with a dual threaded design and ranges in size from 20 to 32mm in length. The RetroFusion screw is inserted between the proximal and middle phalanges, so the opposing threads fixate on the phalangeal canal of the toe and compress the joint.

The TRIM-IT Drill Pin® is a disposable kit containing a bioabsorbable implant, K-Wire, Tamp and Guide Sleeve. The implant has a diameter of 1.5 mm or 2 mm and length of 100 mm. It is pointed at both ends.

The TRIM-IT Spin Pin™ is a disposable kit containing an implant and ruler. The implant has a diameter of 2 mm with a length of 102 mm. One end is made of bioabsorbable material, while the other end has a metal tip. The metal tip has graduated lines for measuring. Both ends are pointed.

The TRIM-IT™ Fixation Screws are fully threaded, cannulated bioabsorbable implants available in various diameters. The unique characteristic of the TRIM-IT screw is the ability of the surgeon to cut the screw to the patient appropriate length.

B. INDICATIONS

The PIP Darts and RetroFusion screws are intended to fix small bony or apical chondral fragments in the foot, ankle, upper extremities, hand, and wrist, where such fragments are not under tension or load-bearing. These devices are used in cases of osteochondritis dissecans and osteochondral fragments, fixation of fractures, 1° metatarsal (bunionectomy osteotomies), cuneiform bones, inherently stable osteotomies, and fusions of the phalanges, metatarsals, metacarpals, carpal bones, tarsal bones, ankle, and wrist. The devices can be used for inherently stable intramedullary stabilization of joint arthroplasty (resection) or fusion for the treatment of digital deformities of the foot or hand. This device is also used in inherently stable long bone fractures such as the femur, fibula, tibia, radius and ulna, including the diaphyseal, epiphyseal, and metaphyseal areas.

The Drill Pin and the Spin Pin are intended for fixation of small bony or apical chondral fragments in the Foot and Ankle where such fragments are not under tension or load-bearing. These pins are used in cases of osteochondritis dissecans and osteochondral fragments, fixation of fractures, 1° metatarsal (bunionectomy osteotomies), cuneiform bones, inherently stable osteotomies, and fusions of the phalanges, metatarsals, tarsal bones, and ankle. The Drill Pin and the Spin Pin can be used for inherently stable intramedullary stabilization of joint arthroplasty (resection) or fusion for the treatment of digital deformities of the foot.

The TRIM-IT Fixation Screw System is a family of implants for the maintenance of alignment and fixation of fractures, osteotomies, arthrodeses or condylar grafts of the foot and ankle in the presence of an appropriate brace and/or immobilization. More specific surgeries include:

Foot/Ankle: Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction (proximal and chevron), Midfoot Reconstruction, Metatarsal fusions and Ligament Repair, Tendon transfer in the foot/ankle such as Posterior Tibial Tendon Transfer for Posterior Tibial Tendon Dysfunction, Flexor Digitorum Longus Transfer for Posterior Tibial Tendon Dysfunction, Kidner transfer, Extensor Hallucis Longus Transfer, Flexor Hallucis Longus for Achilles Tendon Reconstruction, Weber A and B fractures, Lateral and Medial Malleolus ankle Fracture, Osteochondritis Dissecans, Fractures of the Distal Radius, tibial tubercle avulsions in Adolescence

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may tend to retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infection, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

E. WARNINGS

- An internal fixation device must never be reused.
- Only metal implants can be re-sterilized.
- Over-flexing of the device may cause breakage or cracking.
- Postoperatively, until healing is complete the fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the implant, otherwise loosening, fracture, or migration of the device may result.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Metal only:** Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- The 1.5mm Drill Pin can be used with the hammer toe indication only.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

F. MRI SAFETY INFORMATION

I. MR CONDITIONAL

Non-clinical testing and in-vitro electromagnetic simulations demonstrated that the metal (titanium and stainless steel) PIP Darts, RetroFusion Screws and TRIM-IT Devices are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

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