English

A. DEVICE DESCRIPTION

The Jones Bolt Nail and Compression Head are designed to provide fixation of fifth metatarsal fractures while they heal.

B. INDICATIONS

The Jones Bolt Nail Fifth Metatarsal System is intended for use in the fixation of fractures and osteotomies of the fifth metatarsal, including Jones, avulsion and shaft fractures.

Caution: Federal law restricts this device to sale by or on the order of a physician.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone/soft tissue.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation. Nickel sensitivity: If suspected, tests are performed prior to implantation.
- Any active or latent infection or blood supply limitations.
- Patients with poor prognosis for healing.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Potential adverse events related to any surgery: Reactions to anesthesia, the anesthetic or other medications, bleeding infection (both deep and superficial), ileus, blood vessel damage, nerve or soft tissue damage, atelectasis, pneumonia, hematoma, seroma, wound dehiscence or incisional hernia, urologic problems, embolism, anemia, colitis, thrombophlebitis, heart attack, stroke or death.
- Potential adverse events related to the device: Bending, cracking, or fracture of the implant, failure to achieve fracture healing, or loosening of the implant. Limb shortening or loss of anatomic position with nonunion, malunion with rotation or angulation. Irritational injury of soft tissues, including impingement syndrome, or tissue reactions which include macrophage and foreign body reactions adjacent to implants, Metal sensitivity reactions and/or allergic reactions to foreign materials, pain, discomfort, or abnormal sensations due to the presence of the implant. Necrosis of bone or bone resorption. Necrosis of the tissue or inadequate healing. Generation of material debris possibly resulting in bodily response. Additional surgery may be necessary for implant removal, repositioning or replacement.

E. WARNINGS

- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- The correct selection and placement of the implant is extremely important. The appropriate type and size should be selected for the patient. Failure to use the correct implant size or improper positioning may result in loosening, bending, cracking, or fracture of the device, or bone or both. The patient must be cautioned about the use, limitations, and possible adverse effects of this implant including the possibility of the device failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight and/or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device, the patient leaflet and the patient implant card should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Nail and Compression Head:** Removal of supplemental fixation after healing. If the supplemental fixation is not

removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

- Follow your institutions policy for safe disposal of all needles and other sharps or medical waste.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and electromagnetic simulations demonstrated that the Jones Bolt System is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Jones Bolt System is expected to produce a maximum temperature rise less than 6 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Jones Bolt System extends approximately 40 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Implants are single-use only and should never be reused.
- Use care in handling and storage of implant components. Cutting, sharply bending or scratching the surface can significantly reduce the strength and fatigue resistance of the implant system.
- Implants should not be forced, twisted, bent, or pried during insertion as this may damage, break, or bend the implant and render it unusable. Damaged implants should not be implanted.
- The surgeon should take care with patient selection when treating patients with preexisting conditions that may affect the success of the surgical procedure.
- The patient should be advised that an additional procedure for the removal of the implant is sometimes necessary. Providing postoperative instructions to patients and appropriate nursing care is critical. Early weight and/or load bearing substantially increases implant loading and increases the risk of loosening, bending or breaking the device. Patients who are obese and/or noncompliant, as well as patients who could be pre-disposed to delayed union or non-union, must have auxiliary support. Patients should be cautioned against unassisted weight-bearing activity without physician direction or medical release. Postoperative care and physical therapy should be structured to prevent loading of the operative extremity until directed by the physician.
- Even after full healing, the patient should be cautioned that re-fracture is possible and that a revision surgery may be necessary.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Stainless Steel

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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