English

A. DEVICE DESCRIPTION

The Interference Screw family includes the Soft Screw™ , Interference Screw, Tenodesis Screw™ , RetroScrew® , Round Delta Interference Screw, Delta Tapered Interference Screw, FastThread™ Interference Screw and Bio-Cortical™ Screw.

These interference screw products may be threaded or partially threaded, and cannulated or non-cannulated.

B. INDICATIONS

Interference Screws, RetroScrew® and Bio-Cortical™ Screws: These devices are intended to be used for fixation of tissue, including ligament or tendon to bone, or a bone/tendon to bone. See below for specific indications. Tenodesis Screw™: Intended to provide soft tissue reattachment (i.e. fixation of ligament and tendon graft tissue).

Note: Surgeons must apply their professional judgment when determining the appropriate screw size based on the specific indication, preferred surgical technique, and patient history.

See below for specific indications.

Shoulder (Interference Screws only, except Delta Tapered): Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromion-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulolabral Reconstruction.

Foot/Ankle (except Delta Tapered screw): Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Midfoot Reconstruction, Metatarsal Ligament Repair, Flexor Hallucis Longus for Achilles Tendon reconstruction, tendon transfers in the foot and ankle

Knee: Anterior Cruciate Ligament Repair, Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis, **(BioComposite Interference Screw, FastThread™ Interference Screw, Round Delta BioComposite Interference Screw, Delta Tapered BioComposite Interference Screw, BioComposite RetroScrew)** Posterior Cruciate Ligament Repair

Elbow (except Delta Tapered screw): Biceps Tendon Reattach-ment, Ulnar or Radial Collateral Ligament Reconstruction

Hand/Wrist (except Delta Tapered screw): Scapholunate Ligament Reconstruction, Ulnar Collateral Ligament Reconstruction, Radial Collateral Ligament Reconstruction, Carpometacarpal joint arthroplasty (basal thumb joint arthroplasty), Carpal Ligament Reconstructions and repairs, tendon transfer in the hand/wrist

C. CONTRAINDICATIONS

- Bio-Tenodesis only:** Screw sizes that are smaller than 7mm may not be appropriate for the knee indication.
- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Metal implants only:** Shoulder dislocation/subluxation.

E. WARNINGS

- An internal fixation device must never be reused.
- Metal implants only:** All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Bioabsorbable implants:** Do not re-sterilize this device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Metal implants only:** Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Metal Implants Only:** Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

F. MRI SAFETY INFORMATION

I. MR CONDITIONAL

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Interference Screws are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less

- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Interference Screws are expected to produce a maximum temperature rise of up to 14.9° C after 15-minutes of continuous scanning.

- Artifact Information**

In non-clinical testing, the image artifact caused by the Interference Screws can extend up to approximately 120 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

II. MR SAFE

The Interference Screws devices manufactured from only polyetheretherketone (PEEK), polypropylene or poly (L-Lactide), poly (L, DL-Lactide; PLDLA), biphasic calcium phosphate, biphasic tricalcium phosphate (TCP) and/or hydroxyapatite (HA) are MR safe.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Under insertion of the device may leave the proximal end of the implant protruding beyond the cortical bone, which could potentially cause soft tissue irritation and/or pain post-operatively.
- Bio-Tenodesis Screw only:** Use the appropriate size Arthrex Drill to create a pilot hole in the bone.
- Bio-Cortical and Delta Tapered Interference Screw only:** Insert the screwdriver into the screw to a fully seated position. Failure to engage the screw completely may result in damage to the implant.
- Bio-absorbable Interference Screw only:** It is important to completely seat the screwdriver to prevent potential stripping of the hex and/or screw fracture during insertion or removal.
- Bio-absorbable Interference Screw only:** If inserting the Interference Screw through the anteromedial portal, a knee flexion angle of 120° must be maintained throughout the entire insertion process. **Not maintaining or changing the knee flexion angle during screw insertion may result in screw divergence or failure of the screwdriver.** If achieving and maintaining the appropriate flexion angle is not possible or reasonable, a central transpatellar tendon portal should be considered for proper insertion.
- Bio-absorbable Interference Screw only:** During screw size changes or revision procedures, the sheath may be reinserted with the screwdriver over a guide pin to back out the screw in the sheath for arthroscopic removal.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbolsglossary**.

I. STERILIZATION

This device is provided sterile or non-sterile. Refer to the package label for the sterilization method.

Metal implants only: This device can be re-sterilized. It must be adequately cleaned, then sterilized using one of the following sterilization parameters.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
Gravity-Displacement Steam Sterilization Cycle	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C (270°F)	15 Minutes	15 to 30 Minutes
	135°C (275°F)	10 Minutes	30 Minutes
Pre-vacuum Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

STERILIZATION PARAMETERS- FOR OUTSIDE THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
Gravity-Displacement Steam Sterilization Cycle	132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes

Certain Arthrex devices that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI S179, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities” for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions.

Cooling – The device must be adequately cooled after being removed from the sterilizer.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The device is made of titanium, polyetheretherketone (PEEK), polypropylene or poly (L-Lactide). Bio-composite devices are made of made of various configurations of polymers including poly (L-lactide), poly (L, DL-lactide; PLDLA), biphasic calcium phosphate, biphasic tricalcium phosphate (TCP) and hydroxyapatite (HA).

Sheath (if present): Polypropylene, silicone dam

K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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