English

A. DEVICE DESCRIPTION

The Interference Screw family includes the Soft Screw™ , Interference Screw, Tenodesis Screw, RetroScrew®, Round Delta Interference Screw, Delta Tapered Interference Screw, FastThread™ Interference Screw and Bio-Cortical™ Screw.

These interference screw products may be threaded or partially threaded, and cannulated or non-cannulated.

B. INDICATIONS

Interference Screws, RetroScrew, and Bio-Cortical Screws: These devices are intended to be used for fixation of tissue, including ligament or tendon to bone, or a bone/tendon to bone. See below for specific indications.

Tenodesis Screw: Intended to provide soft tissue reattachment (i.e. fixation of ligament and tendon graft tissue).

Note: Surgeons must apply their professional judgment when determining the appropriate screw size based on the specific indication, preferred surgical technique, and patient history.

See below for specific indications.

Shoulder (Interference Screws only, except Delta Tapered): Acrotor Clav Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Rotator-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulolabral Reconstruction

Foot/Ankle (except Delta Tapered screws): Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Midfoot Reconstruction, Metatarsal Ligament Repair, Flexor Hallucis Longus for Achilles Tendon reconstruction, tendon transfers in the foot and ankle

Knee: Anterior Cruciate Ligament Repair *(except 2.5 mm Tenodesis Screw), Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis, **(BioComposite Interference Screw, FastThread Interference Screw, Round Delta BioComposite Interference Screw, Delta Tapered BioComposite Interference Screw, BioComposite RetroScrew)** Posterior Cruciate Ligament Repair*

Elbow (except Delta Tapered screw): Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction

Hand/Wrist (except Delta Tapered screw): Scapholunate Ligament Reconstruction, Ulnar Collateral Ligament Reconstruction, Radial Collateral Ligament Reconstruction, Carpal Ligament Reconstructions and repairs, tendon transfer in the hand/wrist, and Carpometacarpal joint arthroplasty (basal thumb arthroplasty) *(except 2.5 mm Tenodesis Screw)*

C. CONTRAINDICATIONS

- Bio-Tenodesis only:** Screw sizes that are smaller than 7mm may not be appropriate for the knee indication.
- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Metal implants only:** Shoulder dislocation/subluxation.

E. WARNINGS

- An internal fixation device must never be reused.
- Metal implants only:** All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Bioabsorbable implants:** Do not re-sterilize this device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Metal implants only:** Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Metal Implants only:** Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

F. MRI SAFETY INFORMATION

I. MR CONDITIONAL

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Interference Screws are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less

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