# English

In accordance with EN ISO 17664 during the validation process limit values and a means for monitoring chemical residues.

The following cleaning have been established for the product and has been analyzed utilizing cytotoxicity testing as the primary means of assessing safety of the detergent and/or disinfectant levels after the cleaning and washing processes have been carried out.

The quality of water utilized for the cleaning and disinfection process is in the following order:

- ANSI/AMI ST79 - Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities
- ISO 17664: Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices

The instructions in this DFU were developed using the guidance given in the following documents:

- ANSI/AMI ST79 - Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities
- ISO 17664: Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
- ISO 10993-5: Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity

The **Arthrex Low Profile Screws (2.0-3.0 mm solid)** are intended to be used as stand-alone bone screws, or in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, carpalus, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used with a plate, the screws may be used with the Arthrex Low Profile, Small Fragment Plates, Distal Extremity Plates, Mesh Plates and Distal Radius Plates.

The **Arthrex Low Profile Screws (2.0-3.0 mm cannulated)** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, and wrist.

The **Arthrex Low Profile Screws (3.5 mm larger, solid)** are intended to be used as stand-alone bone screws, or in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used with a plate, the screws may be used with the Arthrex Low Profile, Small Fragment Plates, Humeral Fracture Plates, and Osteotomy Plates.

The **Arthrex Low Profile Screws (3.5 mm and larger, cannulated)** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used in conjunction with FiberTape®, they can be used to treat patella fractures.

The **Arthrex Blunt Tip Screws** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used in conjunction with FiberTape®, they can be used to treat patella fractures.

The **Arthrex Vented Cannulated Locking Screws (4.0 mm, cannulated)** are intended to be used in a plate-screw system for internal bone fixation for bone fractures in the humerus. The screws are used with the Arthrex Humeral Fracture Plates. These screws are for **EU distribution only**.

The **Arthrex Blunt Tip Screws** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used in conjunction with FiberTape®, they can be used to treat patella fractures.

The **Arthrex Vented Cannulated Locking Screws (4.0 mm, cannulated)** are intended to be used in a plate-screw system for internal bone fixation for bone fractures in the humerus. The screws are used with the Arthrex Humeral Fracture Plates. These screws are for **EU distribution only**.

The **Arthrex Blunt Tip Screws** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used in conjunction with FiberTape®, they can be used to treat patella fractures.

- Insufficient quantity or quality of bone.
- Foreign body sensitivity.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

#### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

#### E. WARNINGS

- In internal fixation devices must never be reused.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device. Reuse of any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

Phase	Recirculation Time	Temperature	Detergent
<b>Cold Pre-wash</b>	3 Minutes	68 ± 9°F (20 ± 5°C)	N/A
<b>Cleaning Wash</b>	10 Minutes	150 ± 9°F (65.5 ± 5°C)	Enzymatic or alkaline agent
<b>Rinse 1</b>	3 Minutes	122 ± 9°F (50 ± 5°C)	N/A
<b>Rinse 2</b>	3 Minutes	122 ± 9°F (50 ± 5°C)	N/A
<b>Thermal Disinfection Rinse</b>	5 Minutes	194 ± 9°F (90 ± 5°C)	N/A
<b>Drying</b>	6 Minutes	239°F (115°C)	N/A

#### F. MRI SAFETY INFORMATION

- This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for healing, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Use the appropriately sized drill bit for the screw.
- Damage to the driver or screw may result from failure to seat the driver fully into the screw or to align the driver properly with the screw.
- QuickFix™ Screws Only:** Bending the insertion post to remove it from the screw head is not recommended. Screws should be inserted by hand and not with powered equipment.

- PACKAGING AND LABELING**
  - Arthrex devices should be accepted only if the factory packaging and labeling have not been altered.
  - Contact Customer Service if the package has been opened or altered.
  - All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

#### I. VALIDATION

The recommended cleaning, disinfection, and sterilization methods in this DFU have been validated in compliance with federal and international guidance/standards. In accordance with ISO 17665, the half cycle "overkill" approach was used for sterilization validation, and demonstrates a sterility assurance level (SAL) of 10<sup>-6</sup>.

Cleaning, disinfecting, and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility/end user to use equipment operators in the FDA or the local governing body approved cleaning and disinfection with directions from this DFU.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: <i>FOR THE USA ONLY:</i>			
<b>Exposure Temperature</b>	<b>Exposure Time</b>	<b>Drying Time</b>	
121°C (250°F)	30 Minutes	15 to 30 Minutes	
<b>Gravity-Displacement Steam Sterilization Cycle</b>	132°C (270°F)	15 Minutes	15 to 30 Minutes
	135°C (275°F)	10 Minutes	30 Minutes
<b>Pre-vacuum Cycle</b>	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

**STERILIZATION PARAMETERS: *FOR OUTSIDE THE USA ONLY:***

<b>Exposure Temperature</b>	<b>Exposure Time</b>	<b>Drying Time</b>	
132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes	
<b>Gravity-Displacement Steam Sterilization Cycle</b>	121°C (250°F)	30 Minutes	15 to 30 Minutes
<b>Pre-vacuum Cycle</b>	132°C – 135°C (270°F – 275°F)	4 Minutes	20 to 30 Minutes

Cooling – The device must be adequately cooled, after being removed from the sterilizer. It should not be touched during the cooling process. Do not place the device on a cold surface or immerse in a cold fluid.

#### L. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. This device is made of titanium or stainless steel.

#### M. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

#### N. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

- Immerse devices in disinfection solutions for a minimum of 20 minutes. Suitable disinfection solutions can include, but are not limited to: CIDEX®, WAVECIDE®-01, Gaseps®, Kohlsoln®, and equivalent products). Use the supplier's instructions for preparing the solution. **CAUTION: Low acid or high alkaline solutions are not recommended, as these corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™.** If non-neutral pH disinfection chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as not to negatively impact the fit, finish, or function of the device.
- Disinfection solutions should always be mixed to the manufacturer's specification for concentration.
- After disinfection, the devices should be rinsed with cold distilled water or deionized sterile water. (Acceptable an end use of deionized water per AMI TR 34:2014 is <10 EU/ml).
- Dry devices thoroughly utilizing compressed air, wipes, or an oven.
- Check devices for visible soil. Repeat cleaning if soil is visible and re-inspect.

#### II. MACHINE (AUTOMATIC) CLEANING, AND THERMAL DISINFECTION

The instructions for Machine (Automatic) Cleaning and Thermal Disinfection may be followed as an alternative cleaning method to Manual Cleaning and Disinfection.

- Load the devices in the washer such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain.
- Run the automatic wash cycle. The minimum cycle parameters are listed below:

MINIMUM WASHING CYCLE PARAMETERS			
Phase	Recirculation Time	Temperature	Detergent
<b>Cold Pre-wash</b>	3 Minutes	68 ± 9°F (20 ± 5°C)	N/A
<b>Cleaning Wash</b>	10 Minutes	150 ± 9°F (65.5 ± 5°C)	Enzymatic or alkaline agent
<b>Rinse 1</b>	3 Minutes	122 ± 9°F (50 ± 5°C)	N/A
<b>Rinse 2</b>	3 Minutes	122 ± 9°F (50 ± 5°C)	N/A
<b>Thermal Disinfection Rinse</b>	5 Minutes	194 ± 9°F (90 ± 5°C)	N/A
<b>Drying</b>	6 Minutes	239°F (115°C)	N/A

- Automatic wash cleaning solutions can include, but are not limited to: neidisher® Mediclean forte, Thermosept® alkal clean, Prolystic® Ultra Concentrate Enzymatic Cleaner, and Proklenz NpI Neutral Detergent. **CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylene-propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™.** If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.
- Check devices for visible soil. Repeat cleaning if soil is visible and re-inspect.

#### K. STERILIZATION

This device is provided sterile or non-sterile. Check the package labeling for more information. For the devices that are not provided in a terminally sterilized configuration, sterilization is to be performed following cleaning, disinfection prior to use.. This device (i.e. for open but unused) may be re-sterilized following cleaning, disinfection prior to use.

Sterilization may be performed using the specified sterilization cycles shown below. Ensure that the wrap is large enough to contain the device without stressing the seals.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to the Arthrex Instruments Directions for Use and ANSI/AMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.

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